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IMPROVE REVIEW

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Thromboelastography

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

- 赵绥民,田兆嵩.急性失血病人的成分输血[J].《中国输血杂志》,1999;12(1):45-46.
- 单闯.失血性休克的病理生理及治疗[J].《现代实用医学》,2002;14(11):620-622.
- KM Comer. Anemia as a feature of primary gastrointestinal neoplasia[J]. Compendium on Continuing Education for the Practicing Veterinarian Immun. Microbiol. infect. Dis. 1990; 13(4): 217-221.
- Norton R., Kobusingye O. Injuries[J]. N Engl J Med. 2013; 368(18): 1723-1730.
- Rhee P., Joseph B., Pandit V., Aziz H., Vercruyse G., Kulvatunyou N., Friese R. S. Increasing trauma deaths in the United States[J]. Ann Surg. 2014; 260 (1): 13-21.
- PI Johansson, AM Sørensen, CF Larsen, NA Windeløv, J Stensballe, A Perner, LS Rasmussen, SR Ostrowski. Low hemorrhage-related mortality in trauma patients in a Level I trauma center employing transfusion packages and early thromboelastography-directed hemostatic resuscitation with plasma and platelets[J]. Transfusion. 2013; 53(12):3088-3099.
- Holcomb JB, Jenkins D, Rhee P, et al. Damage control resuscitation: directly addressing the early coagulopathy of trauma[J]. Trauma. 2007; 62(2): 307-310.
- RM Durham, JJ Moran, JE Mazuski, MJ Shapiro, AE Baue, LM Flint. Multiple organ failure in trauma patients[J]. The Journal of trauma. 2003; 55(4): 608-616.
- Cohen MJ. New understandings of post injury coagulation and resuscitation[J]. International journal of surgery (London, England). 2016; 33:242-245.
- Satoshi Gando, Yasuhiro Otomo. Local hemostasis, immunothrombosis, and systemic disseminated intravascular coagulation in trauma and traumatic shock [J]. Critical Care. 2015; 19: 72.
- Gando S. Pathophysiology of Trauma-Induced Coagulopathy and Management of Critical Bleeding Requiring Massive Transfusion[J]. Semin Thromb Hemost. 2016; 42: 155-165.
- L.G. Stansbury, A.S. Hess, K. Thompson, et al. The clinical significance of platelet counts in the first 24 hours after severe injury[J]. Transfusion. 2013; 5:34.
- Kutcher ME, Redick BJ, McCreery RC, et al. Characterization of platelet dysfunction after trauma[J]. The journal of trauma and acute care surgery, 2012; 73(1):13-19.
- Sillesen M, Johansson PI, Rasmussen LS, et al. Platelet activation and dysfunction in a large-animal model of traumatic brain injury and hemorrhage[J]. The journal of trauma and acute care surgery. 2013; 74(5): 1252-1259.
- M.V. Wohlauer, E.E. Moore, S. Thomas, et al. Early platelet dysfunction: an unrecognized role in the acute coagulopathy of trauma[J]. J. Am. Coll. Surg. 2012; 214: 739-746.
- 刘卓刚,胡荣.我国不同人群缺铁性贫血的流行病学特点[J].中国实用乡村医生杂志,2008;15(6):5-6.
- 中华医学会.临床诊疗指南·血液学分册[M].北京:人民卫生出版社,2009:1-3.
- 黄才千,仇小强.儿童缺铁性贫血预防策略研究[J].中华疾病控制杂志,2008;12(4):391-394.
- Meng XD, Xiang GQ, He XH, et al. Etiology and clinical character of liver cirrhosis: an analysis of cases[J]. World J Digestol, 2008; 16(7): 1880-1884.
- Rodriguez-Castro KI, Antonello A, Ferrarese A. Spontaneous bleeding or thrombosis in cirrhosis: What should be feared the most?[J]. World J Hepatol, 2015; 7(14): 1818-1827.
- Tripodi A, Mannucci PM. The coagulopathy of chronic liver disease[J]. N Engl J Med, 2011; 365(2): 147-156.
- Cho J, Choi SM, Yu SJ, et al. Bleeding complications in critically ill patients with liver cirrhosis[J]. Korean J Intern Med, 2016; 31(2): 288-295.
- Andriulli A, Tripodi A, Angeli P, et al. Hemostatic balance in patients with liver cirrhosis: Report of a consensus conference [J]. Dig Liver Dis, 2016; 48(5): 455-467.
- Jairath V, Burroughs AK. Anticoagulation in patients with liver cirrhosis: complication or therapeutic opportunity? [J]. Gut, 2013; 62(4): 479-482.
- Marten M, Coll M, Ezkurdia N, et al. Physiopathology of splanchnic vasodilation in portal hypertension[J]. World J Hepatol, 2010; 2(6): 208-220.
- Garbuzenko DV. The mechanisms of adaptation of the vascular bed to hemodynamic changes in portal hypertension [J]. Vestn Ross Akad Med Nauk, 2013; (1): 52-57.

27. 卢尧, 滕方, 孙桂香等. 肝硬化患者凝血障碍机制研究进展[J]. 中国输血杂志, 2015, 28(2): 216-220.
28. Iwakiri. Pathophysiology of portal hypertension[J]. Clin Liver Dis, 2014; 18(2): 281-291.
29. 减立娜, 刑慧慧, 王嘉菲等. 肝硬化并发门静脉血栓研究进展[J]. 实用肝脏病学杂志, 2014; 17(4): 429-432.
30. Roberts LN, Bernal W. Management of Bleeding and Thrombosis in Critically Ill Patients with Liver Disease[J]. Semin Thromb Hemost, 2015; 41(5): 520-526.
31. 胡晓, 冯立民, 王佳丽. 肝硬化患者血浆凝血因子VIII及常规凝血指标检测结果分析[J]. 实验与检验医学, 2013; 31(1): 24-26.
32. Tapper EB, Robson SC, Malik R. Coagulopathy in cirrhosis the role of the platelet in hemostasis[J]. J Hepatol, 2013; 59(4): 889-890.
33. Gand R, Chunlei H. Study on the relationship between coagulation indexes in patients with liver cirrhosis[J]. Int J Lab Med, 2013; 24(32): 3019-3020.
34. Zongbao Liu, Erqing Chai, Hecheng Chen, et al. Comparison of Thrombelastography (TEG) in Patients with Acute Cerebral Hemorrhage and Cerebral Infarction[J]. Medical Science Monitor. 2018; 24: 6466-6471.
35. Zahr Eldeen F, Roll GR, Derosas C, et al. Preoperative thrombelastography as a sensitive tool predicting those at risk of developing early hepatic artery thrombosis after adult liver transplantation[J]. Transplantation. 2016; 100(11): 2382-2390.
36. Chitlur M, Warriar I, Rajpurkar M, et al. Thrombelastography in children with coagulation factor deficiencies[J]. Br J Haematol, 2008; 142 (2) : 250-256.
37. Zia AN, Chidur M, Rsipurkar M, et al. Thrombelastography identifies children with rare bleeding disorders and predicts phenotype[J]. Haemophilia, 2015; 21(1) : 124-132.
38. Müller MC, Meijers JC, Vroom MB, et al. Utility of thrombelastography and/or thrombelastometry in adults with sepsis: a systematic review[J]. Crit Care, 2014; 18(1): R30.
39. 孙云雨, 朱永亮, 华岚等. 血栓弹力图指导主动脉夹层患者术后成分输血作用的探讨[J]. 临床血液学杂志, 2015; 28(12): 1022-1024.
40. Girdauskas E, Kempfert J, Kuntze T, et al. Thrombelastometrically guided transfusion protocol during aortic surgery with circulatory arrest: a prospective, randomized trial[J]. J Thorac Cardiovasc Surg, 2010; 140: 1117-1124.
41. Ferraris VA, Saha SP, Oestreich JH, et al. 2012 update to the Society of Thoracic Surgeons guideline on use of antiplatelet drugs in patients having Cardiac and noncardiac operations[J]. Ann Thoracic Surg, 2012; 94(5): 1761-1781.
42. Krzanicki D, Sugavanam A, Mallets S. Intraoperative hypercoagulability during liver transplantation as demonstrated by thrombelastography[J]. Liver Transpl, 2013; 19(8): 852-861.
43. Pekelharing J, Furck A, Banya W, et al. Comparison between thrombelastography and conventional coagulation tests after cardiopulmonary bypass surgery in the paediatric intensive care unit [J]. Int J Lab Hematol, 2014; 36 (4) ; 465-471.
44. Bischof D, Dalbert S, Zollinger A, et al. Thrombelastography in the surgical patient[J]. Minerva Anesthesiol, 2010; 76: 383.
45. Butwick A, Ting V, Ralls LA, et al. The association between thrombelastographic parameters and total estimated blood loss in patients undergoing elective cesarean delivery[J]. Anesth Analg, 2011; 112(5): 1041-1047.
46. AJ Wikkelsøe, A Afshari, J Stensballe, et al. The FIB-PPH trial: fibrinogen concentrate as initial treatment for postpartum haemorrhage: study protocol for a randomised controlled trial[J]. Trials, 2012; 13(1): 1-8.
47. Karlsson O, Jeppsson A, Hellgren M. Major obstetric haemorrhage: monitoring with thrombelastography, laboratory analyses or both[J]. Int J Obstet Anesth, 2014; 23(1): 10-17.
48. Welsh KJ, Padilla A, Dasgupta A, et al. Thrombelastography is a suboptimal test for determination of the underlying cause of bleeding associated with cardiopulmonary bypass and may not predict a hypercoagulable state[J]. Am J Clin Pathol, 2014; 142(4): 492-497.
49. 赵凌, 张敏, 徐智杰, 等. 肝素酶纠正的血栓弹力图在胸心外科患者术后的应用[J]. 临床检验杂志, 2014; 32(11): 825-827.
50. 燕姿辰, 杨新春, 张新等. 用血栓弹力图评价肝素浓度变化对凝血系统的影响[J]. 沈阳药科大学学报, 2012; 29(8): 644-649.
51. Tantry US, Bonello L, Aradi D, et al. Consensus and update on the definition of on treatment platelet reactivity to adenosine diphosphate associated with ischemia and bleeding[J]. J Am Coll Cardiol, 2013; 62(24): 2261-2273.
52. 马璐璐, 虞雪融, 黄宇光. 快速血栓弹力图监测关节置换术患者围术期的凝血功能[J]. 中国输血杂志, 2015; 28(6): 655-657.
53. Coppell JA, Thalheimer U, Zamboni A, et al. The effect of unfractionated heparin, low molecular weight heparin and danaparoid on the thrombelastogram (TEG): an in-vitro comparison of standard and heparinase-modified TEGs with Conventional coagulation assays[J]. Blood coagulat fibrinol, 2006; 17(2): 97-104.
54. 纪宏文, 马丽, 高旭蓉等. 血栓弹力图在体外循环心血管手术的应用[J]. 中国体外循环杂志, 2011, 9(3): 170-172.
55. Brohi K, Singh J, Heron M, Coats T. Acute traumatic coagulopathy[J]. The Journal of trauma 2003; 54: 1127-1130.
56. MacLeod JB, Lynn M, McKenney MG, Cohn SM, Martha M. Early coagulopathy predicts mortality in trauma[J]. The Journal of trauma. 2003; 55: 39-44.
57. Maegele M, Lefering R, Yucel N, et al. Early coagulopathy in multiple injury: an analysis from the German Trauma Registry on 8724 patients[J]. Injury 2007; 38: 298-304.
58. Brohi K, Cohen MJ, Ganter MT, Matthay MA, Mackersie RC, Pittet JF. Acute traumatic coagulopathy: initiated by hypoperfusion: modulated through the protein C pathway?[J]. Annals of surgery 2007; 245: 812-818.
59. Niles SE, McLaughlin DF, Perkins JG, et al. Increased mortality associated with the early coagulopathy of trauma in combat casualties[J]. The Journal of trauma. 2008; 64: 1459-63; discussion 63-65.

60. Kaufmann CR, Dwyer KM, Crews JD, Dols SJ, Trask AL. Usefulness of thrombelastography in assessment of trauma patient coagulation[J]. The Journal of trauma. 1997; 42:716-20; discussion 20-22.
61. Park MS, Martini WZ, Dubick MA, et al. Thrombelastography as a better indicator of hypercoagulable state after injury than prothrombin time or activated partial thromboplastin time[J]. The Journal of trauma. 2009; 67:266-75; discussion 75-76.
62. Müller MC, Meijers JC, Vroom MB, et al. Utility of thrombelastography and/or thrombelastometry in adults with sepsis: a systematic review[J]. Crit Care, 2014; 18(1):R30.
63. 文爱清, 蒋建新. “严重创伤输血专家共识”解读[J]. 中华创伤杂志, 2013; 29(8): 706-710.
64. Gary JL, Schneider PS, Galpin M, et al. Can Thrombelastography Predict Venous Thromboembolic Events in Patients with Severe Extremity Trauma ?[J]. J Orthop Trauma, 2016; 30(6): 294-298.
65. Cotton BA, Minei KM, Radwan ZA, et al. Admission rapid thrombelastography predicts development of pulmonary embolism in trauma patients[J]. Trauma Acute Care Surg, 2012; 72(6):1470-1477.
66. Carroll RC, Craft RM, Langdon RJ, et al. Early evaluation of acute traumatic coagulopathy by thrombelastography[J]. Translational research:the journal of laboratory and clinical medicine 2009; 154: 34-39.

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1. **PI Johansson, AM Sørensen, CF Larsen, NA Windeløv, J Stensballe, A Perner, LS Rasmussen, SR Ostrowski. Low hemorrhage-related mortality in trauma patients in a Level I trauma center employing transfusion packages and early thromboelastography-directed hemostatic resuscitation with plasma and platelets[J]. Transfusion. 2013; 53(12):3088-3099.**

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
2. **M Maegele, H Schöchl, T Menovsky, H Maréchal, N Marklund, A Buki, S Stanworth. Coagulopathy and haemorrhagic progression in traumatic brain injury: advances in mechanisms, diagnosis, and management[J]. The Lancet. Neurology. 2017; 16(8): 630-647.**

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3. Cohen MJ. New understandings of post injury coagulation and resuscitation[J]. *International journal of surgery (London, England)*. 2016; 33:242-245.

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4. Satoshi Gando, Yasuhiro Otomo. Local hemostasis, immunothrombosis, and systemic disseminated intravascular coagulation in trauma and traumatic shock [J]. *Critical Care*. 2015; 19: 72.

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5. Gando S. Pathophysiology of Trauma-Induced Coagulopathy and Management of Critical Bleeding Requiring Massive Transfusion[J]. *Semin Thromb Hemost*. 2016; 42: 155-165.

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6. L.G. Stansbury, A.S. Hess, K. Thompson, et al. The clinical significance of platelet counts in the first 24 hours after severe injury. *Transfusion*. 2013; 5:34.

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7. Kutcher ME, Redick BJ, McCreery RC, et al. Characterization of platelet dysfunction after trauma [J]. The journal of trauma and acute care surgery, 2012, 73(1):13-19.

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
8. Sillesen M, Johansson PI, Rasmussen LS, et al. Platelet activation and dysfunction in a large-animal model of traumatic brain injury and hemorrhage[J]. The journal of trauma and acute care surgery. 2013; 74(5): 1252-1259.

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10. Meng XD, xiang GQ, He XH, et al. Etiology and clinical character of liver cirrhosis: an analysis of cases [J]. World J Digestol, 2008, 16(7): 1880-1884.

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12. **Cho J, Choi SM, Yu SJ, et al. Bleeding complications in critically ill patients with liver cirrhosis [J]. Korean J Intern Med, 2016, 31(2): 288-295.**

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13. **Andriullia A, Tripodib A, Angeli P, et al. Hemostatic balance in patients with liver cirrhosis: Report of a consensus conference [J]. Dig Liver Dis, 2016, 48(5): 455-467.**

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17. Roberts LN, Bernal W. Management of Bleeding and Thrombosis in Critically Ill Patients with Liver Disease [J]. Semin Thromb Hemost, 2015, 41(5): 520-526.

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18. Gand R, Chunlei H. Study on the relationship between coagulation indexes in patients with liver cirrhosis [J]. Int J Lab Med, 2013, 24(32): 3019-3020.

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19. Zongbao Liu, Erqing Chai, Hecheng Chen, et al. Comparison of Thrombelastography (TEG) in Patients with Acute Cerebral Hemorrhage and Cerebral Infarction [J]. Medical Science Monitor. 2018; 24: 6466-6471.

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20. **Zahr Eldeen F, Roll GR, Derosas C, et al. Preoperative thromboelastography as a sensitive tool predicting those at risk of developing early hepatic artery thrombosis after adult liver transplantation[J]. Transplantation. 2016, 100(11): 2382-2390.**

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21. **Chitlur M, Warriar I, Rajpurkar M, et al. Thromboelastography in children with coagulation factor deficiencies[J]. Br J Haematol, 2008, 142 (2) : 250-256.**

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22. **Zia AN, Chidur M, Rsjpurkar M, et al. Thromboelastography identifies children with rare bleeding disorders and predicts phenotype[J]. Haemophilia, 2015, 21(1) : 124-132.**

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23. Müller MC, Meijers JC, Vroom MB, et al. Utility of thromboelastography and/or thromboelastometry in adults with sepsis: a systematic review[J]. Crit Care, 2014,18(1):R30.

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
24. Girdauskas E, Kempfert J, Kuntze T, et al. Thromboelastometrically guided transfusion protocol during aortic surgery with circulatory arrest: a prospective, randomized trial[J]. J Thorac Cardiovasc Surg, 2010, 140: 1117-1124.

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25. Krzanicki D, Sugavanam A, Mallets S. Intraoperative hypercoagulability during liver transplantation as demonstrated by thromboelastography[J]. Liver Transpl, 2013, 19(8): 852-861.

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26. Pekelharing J, Furck A, Banya W, et al. Comparison between thromboelastography and conventional coagulation tests after cardiopulmonary bypass surgery in the paediatric intensive care unit [J]. Int J Lab Hematol, 2014, 36 (4): 465-471.

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28. Butwick A, Ting V, Ralls LA, et al. The association between thromboelastographic parameters and total estimated blood loss in patients undergoing elective cesarean delivery[J]. *Anesth Analg*, 2011,112(5): 1041-1047.

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29. AJ Wikkelsoe, A Afshari, J Stensballe, et al. The FIB-PPH trial: fibrinogen concentrate as initial treatment for postpartum haemorrhage: study protocol for a randomised controlled trial [J]. *Trials*, 2012,13 (1):1-8.

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30. Karlsson O, Jeppsson A, Hellgren M. Major obstetric haemorrhage: monitoring with thromboelastography, laboratory analyses or both[J]. *Int J Obstet Anesth*, 2014,23(1):10-17.

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31. Welsh KJ, Padilla A, Dasgupta A, et

a1. Thromboelastography is a suboptimal test for determination of the underlying cause of bleeding associated with cardiopulmonary bypass and may not predict a hypercoagulable state[J]. Am J Clin Pathol, 2014, 142(4): 492-497.

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32. Tantry US, Bonello L, Aradi D, et al. Consensus and update on the definition of on treatment platelet reactivity to adenosine diphosphate associated with ischemia and bleeding[J]. J Am Coll Cardiol, 2013, 62(24): 2261-2273.

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34. Gary JL, Schneider PS, Galpin M, et al. Can Thrombelastography Predict Venous Thromboembolic Events in Patients with Severe Extremity Trauma ?[J]. J Orthop Trauma, 2016, 30(6): 294-298.

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35. Cotton BA, Minei KM, Radwan ZA, et al. Admission rapid thrombelastography predicts development of pulmonary embolism in trauma patients[J]. Trauma Acute Care Surg, 2012, 72 (6): 1470-1477.

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
Shiv Kumar Sarin, M.D., D.M., Department of Hepatology, Institute of Liver and Biliary Sciences, D1 Vasant Kunj, New Delhi 110070, India. E-mail: sksarin@ilbs.in.

Current Concepts in Coagulation Profile in Cirrhosis and Acute-on-Chronic Liver Failure

Madhumita Premkumar, M.D., D.M., ^{ID} and Shiv Kumar Sarin, M.D., D.M.

BACKGROUND

Patients with cirrhosis have profoundly altered hemodynamics and hemostatic pathways with procoagulant and anticoagulant mechanisms, resulting in a tenuous “re-balanced” state in the setting of portal hypertension. This balance could be tipped toward either a procoagulant or anticoagulant phenotype by superimposed conditions.¹

Specific coagulation factor (F) V (FV), FVII, FIX, FX, FXI, prothrombin, protein C, and protein S are reduced with concomitant increased FVIII and von Willebrand factor

(vWF) activity in cirrhosis.² Thrombocytopenia, increased nitric oxide, and prostacyclin inhibit platelet function (PF), and higher vWF and FVIII activity support platelet aggregation.³ Thrombocytopenia is the result of splenic sequestration in portal hypertension, decreased hepatic thrombopoietin synthesis, and immune-mediated platelet destruction due to glycoprotein IIb/IIIa platelet surface antigen-antibody interaction triggered by inflammation or sepsis.⁴ Patients with acute liver failure (ALF) have a markedly prolonged international normalized ratio (INR) but preserved thrombin generation potential, and they

Abbreviations: ACLF, acute-on-chronic liver failure; ACT, activated clotting time; AD, acute decompensation; ALF, acute liver failure; ALI, acute liver failure; aPTT, activated partial thromboplastin time; AT, antithrombin; CR, clot rate; CYP P450, cytochrome P450; DIC, disseminated intravascular coagulation; DOAC, direct oral anticoagulant; EACA, ϵ -aminocaproic acid; F, factor; FFP, fresh frozen plasma; IL, interleukin; INR, international normalized ratio; IV, intravenous; LMWH, low-molecular-weight heparin; LT, liver transplantation; Lysis30, lysis index after 30 minutes; MA, maximum amplitude; MELD, Model for End-Stage Liver Disease; PAI, plasminogen activator inhibitor; PC, platelet concentrate; PF, platelet function; PLT, platelet count; POC, point of care; PPC, prothrombin complex concentrate; PRBC, packed red blood cells; PT, prothrombin time; PVT, portal vein thrombosis; rF, recombinant factor; ROTEM, rotational thromboelastometry; SCT, standard coagulation test; SIRS, systemic inflammatory response syndrome; SMVT, superior mesenteric vein thrombosis; SQ, subcutaneous; TEG, thromboelastography; TIPS, transjugular intrahepatic portosystemic shunt; TNF, tumor necrosis factor; tPA, tissue plasminogen activator; TRALI, transfusion-related acute lung injury; ULN, upper limit of normal; VE, viscoelastic; VKA, vitamin K antagonist; vWF, von Willebrand factor.

From the Department of Hepatology, Institute of Liver and Biliary Sciences, New Delhi, India.

Potential conflict of interest: Nothing to report.

Received December 24, 2019; accepted April 22, 2020.

View this article online at wileyonlinelibrary.com

© 2020 by the American Association for the Study of Liver Diseases

tend not to bleed. In cirrhosis, INR increases modestly with synthetic dysfunction and reflects short-term patient mortality, but not bleeding risk, during invasive procedures because of adequate thrombin generation potential.⁵ Acute-on-chronic liver failure (ACLF) is associated with organ failure, progressive clinical course, and high 28-day mortality.⁶ In ACLF, coagulation switches from a procoagulant to an anticoagulant phenotype with the onset of systemic inflammation and endothelial activation, which is compounded by sepsis.⁷⁻⁹ In these situations, critical scenarios, such as variceal bleeding, the need for invasive procedures, central line placement, venous thromboembolism, portal vein thrombosis (PVT), and secondary organ failure, often supervene. Standard coagulation tests (SCTs), such as prothrombin time (PT), activated partial thromboplastin time (aPTT), INR, and bleeding time, do not measure the risk for bleeding in liver disease.³ Hence inappropriate use of blood products to correct these parameters is often ineffective and leads to volume overload and transfusion-related acute lung injury (TRALI).^{10,11}

CELL-BASED MODEL OF HEMOSTASIS

The traditional chemical enzymatic pathways are now superseded by the cell-based model of hemostasis in which activated platelets are the primary effectors of clot initiation making a platelet plug. A crucial cellular element is endothelial injury, which is a driver for bleeding and thrombosis with activation of the inflammatory cascade. Primary hemostasis is mediated by platelets, whereas secondary hemostasis is characterized by formation of the fibrin mesh with plasma procoagulant proteins.⁴ Finally, when the vascular repair is complete, plasma anticoagulant proteins break down this fibrin mesh in a third process known as *fibrinolysis*. In cirrhosis, bleeding is often due to portal hypertension per se rather than being coagulopathy related, like hyperfibrinolysis. The cell-based model also explains why local hemostatic changes at the site of injury do not override the systemic hemostatic balance, and conventional tests of coagulation remain the same in a patient with liver disease with clinically apparent bleeding.^{1,2} The clinician's ability to detect the coagulation defect at the site of portal hypertensive bleeding, such as variceal or mucosal or vessel injury (e.g., after biopsy), remains elusive because conventional tests reveal the global rather than the local defect.¹²

PATHOPHYSIOLOGY OF COAGULATION DEFECTS IN CIRRHOSIS

Patients with compensated liver cirrhosis have a rebalanced coagulation profile. In patients with cirrhosis, bleeding (mostly variceal) is usually the consequence of portal hypertension (gastropathy/variceal bleeding) or due to vessel injury (biopsy site/paracentesis site) bleeding.³ This implied that an elevated INR is not a predictor of bleeding, nor is a therapeutic correction target in a bleeding patient.⁵ Progression of chronic liver failure is associated with derangement of coagulative balance, and the balance shifts from a procoagulant to anticoagulant phenotype. Platelet count (PLT) alone is an incomplete test of coagulation and indirectly correlates with the degree of portal hypertension, splenomegaly, and hepatic decompensation. However, studies have shown that a target PLT of $56 \times 10^9/L$ suffices to control variceal bleeding because of intact platelet-dependent thrombin generation in cirrhosis.^{13,14} A special scenario is liver transplantation (LT) in cirrhosis, where excess transfusion of platelet concentrate (PC), cryoprecipitate, or fresh frozen plasma (FFP) may contribute to the development of hepatic arterial or venous thrombosis.^{15,16}

PATHOPHYSIOLOGY OF COAGULATION IN ACLF

Recent studies have sparked interest in coagulation in ACLF.^{3,7-9,11} Fisher et al.⁸ showed lower thrombin generation in patients with ACLF compared with those with acute decompensation (AD). Bedreli et al.⁹ demonstrated that rotational thromboelastometry (ROTEM) may reduce substitution of coagulation factors in patients with ACLF. Blasi et al.⁷ described the ROTEM profile in ACLF as compared with cirrhosis. Finally, our group prospectively studied the dynamic global coagulation dysfunction in patients with ACLF in relation to the development of systemic inflammatory response syndrome (SIRS) and sepsis, and found that worsening coagulation failure often preceded the development of sepsis and mortality in these critically ill patients.^{3,17} Table 1 summarizes available studies on coagulation in cirrhosis and liver failure.^{3,5,7,8,18-24} Figure 1 shows our understanding of this dynamic model of coagulation in cirrhosis and ACLF.

RISK ASSESSMENT FOR PROCEDURES IN LIVER DISEASE

Preprocedure correction of SCTs is a double-edged sword. Assessment of baseline bleeding risk, limitations of

TABLE 1. ROLE OF VE TESTING IN CIRRHOSIS, ALF, AND ACLF

Author (Year)	No. of Patients	Clinical Condition	Findings	Specific Coagulation Parameters	Remarks
Cirrhosis					
Papatheodoridis et al. (1999) ²²	Cirrhosis: n = 84 30 subjects with sepsis	Cirrhosis	All TEG parameters worsened with sepsis and recovered with resolution	Successful treatment of infection usually restores hemostasis parameters to pre-infection levels in 5 days	Association of coagulation with sepsis
Lloyd-Donald et al. (2017) ²³	Cirrhosis: n = 74	Acutely ill patients with cirrhosis	Delayed clot formation with longer reaction time ($P < 0.01$), longer kinetic time ($P < 0.01$), more acute α angle (65 versus 72.2 degrees; $P < 0.01$), and longer ACT ($P < 0.01$)	Not done	AD of cirrhosis made the TEG hypocoagulable
Bedreli et al. (2016) ⁹	Liver cirrhosis: n = 37	Cirrhosis	ROTEM-guided transfusion before minimally invasive procedures	VE test-based blood transfusion protocols	Supplementation of coagulation factors (PCCs, FFP, and fibrinogen)
LT					
Fayed et al. (2015) ²⁴	n = 100	Decompensated cirrhosis during LT	POC test ROTEM	Thromboelastometry-based regression equation accounted for 63% of PRBC, 83% of FFP, 61% of cryoprecipitate, and 44% of platelet transfusion requirements	POCs have ability to predict transfusion requirements
De Pietri et al. (2016) ¹⁹	n = 261	Cirrhosis during LT	Greater clot strength was associated with a MELD score < 20 , hepatitis C virus, and cholestatic-related cirrhosis		
ALF/ALI					
Agarwal et al. (2012) ⁵	ALF: n = 20	ALF	Significant elevation in plasma levels of FVIII and on vWF and microparticles	TEG tracings were consistent with a hypocoagulable state in 20%, normal in 45%, and hypercoagulable in 35%	Coagulation disturbance in patients with ALF is complex
Stravitz et al. (2012) ²⁰	n = 51	ALF/ALI	Markedly decreased procoagulant FV and FVIII levels	Despite a mean INR of 3.4 ± 1.7 (range 1.5-9.6), mean TEG parameters were normal, and MA was higher in patients with ALF than ALI and correlated directly with venous ammonia concentrations	Five individual TEG parameters were normal in 32 patients (63%); low MA was confined to patients with PLTs $< 126 \times 10^9/L$
ACLF					
Blasi et al. (2018) ⁷	ACLF: n = 36 Decompensated cirrhosis: n = 24	ACLF Grade 1: 13 Grade 2: 13 Grade 3: 10	More hypocoagulable characteristics compared with cirrhosis with AD	Prolonged time to initial fibrin formation and clot formation time, and decreased maximum clot firmness and α angle values	VE testing in cirrhosis and ACLF
Fisher et al. (2018) ⁸	Stable cirrhosis: n = 8 Decompensated cirrhosis: n = 44 ACLF: n = 17	ACLF subjects Grade 1: 2 Grade 2: 9 Grade 3: 6 ACLF, followed up prospectively	Thrombin generation peak and velocity were also significantly lower in ACLF compared with healthy controls	Increased vWF and lower ADAMTS-13 in ACLF	Thrombin generation in ACLF reported Small sample size
Premkumar et al. (2019) ³	ACLF: n = 114 Cirrhosis: n = 25		TEG parameters worsened with SIRS and sepsis	Dynamic coagulation changes predicted sepsis and mortality	Prospective data on individual coagulation factors with 28-day mortality assessment
Goyal et al. (2018) ²¹	ACLF: n = 68 Cirrhosis with AD: n = 53	ACLF	Reduced MA and lysis in patients with ACLF	TEG parameters in patients with ACLF are essentially normal, with the exception of reduced MA	Single time assessment with no assessment of changes with SIRS and sepsis

Abbreviations: ACLF, acute-on-chronic liver failure; AD, acute decompensation; ADAMTS-13, a disintegrin and metalloproteinase with a thrombospondin type 1 motif, member 13; ALF, acute liver failure; ALI, acute liver injury; AT, antithrombin; F, factor; FFP, fresh frozen plasma; INR, international normalized ratio; MELD, Model for End-Stage Liver Disease; PCCs, prothrombin complex concentrates; MA, maximum amplitude; POC, point of care; PT, prothrombin time; ROTEM, rotational thromboelastometry; SIRS, systemic inflammatory response syndrome; TEG, thromboelastography; vWF, von Willebrand factor.


FIG 1 Dynamic coagulation profile changes in cirrhosis and ACLF.

TABLE 2. CLASSIFICATION OF RISK FOR BLEEDING DURING PROCEDURES

High-Risk Procedures	Intermediate-Risk Procedures	Lower-Risk Procedures
Intrabdominal/orthopedic/cardiac surgery	Invasive endoscopy (biliary sphincterotomy, endoscopic ultrasound-guided biopsies, cystogastrostomy)	Paracentesis
Intracranial catheter insertion	TIPS	Thoracentesis
Invasive endoscopy (polypectomy)	Lumbar puncture	Dental extraction
Endoscopic submucosal resection	Transarterial chemoembolization or percutaneous locoregional therapies for hepatocellular carcinoma	Endoscopy (e.g., diagnostic, variceal ligation)
	Percutaneous biopsy of extrahepatic organs and lesions	Cardiac catheterization
	Transjugular liver biopsy	Central line placement
		Hepatic venous pressure gradient measurement

Abbreviation: TIPS, transjugular intrahepatic porto-systemic shunt.

prediction by conventional tests, and the benefits of prophylactic correction are to be weighed (Table 2). Correction of the PT/INR may lead to paradoxical bleeding by increasing blood flow in the collateral beds.²⁵ Standard doses of blood components, such as FFP, can only partially correct SCTs. Common dose ranges are from 15 to 30 mL/kg; however, these data have been extrapolated from critically ill patients with trauma and coagulopathy.^{26,27} Increased volumes of FFP infusion predispose to transfusion-associated circulatory overload and TRALI. For every 100-mL rapid expansion of blood volume, portal pressure increases

by 1 mm Hg, which can subsequently cause portal collateral-related bleeding.²⁸

Prothrombin complex concentrates (PCCs) and recombinant FVIIa (rFVIIa) have been designed for specific situations. The use of cryoprecipitate, a small-volume product having vWF, fibrinogen, and fibronectin without need for cross-matching, is an attractive option.

PCCs are available as three-factor (II, IX, X) and four-factor (II, IX, X, with the addition of VII) products. These products have additional factors (protein C, protein S,

antithrombin [AT] III) with or without heparin. The use of PCCs has been standardized only for vitamin K antagonist (VKA)-treated patients, and the dose is based on body weight, INR, and FIX content. rFVIIa has been studied in acute variceal bleeding; however, no mortality benefit was demonstrated²⁹ (Table 3). The thrombopoietin receptor agonist eltrombopag can increase PLTs to a variable degree after a 10- to 20-day course but bears a risk for PVT.³⁰ Two new congeners, avatrombopag and lusutrombopag, are reportedly bereft of this complication, making them of clinical interest.^{31,32}

EVIDENCE-BASED ROLE OF VISCOELASTIC TESTING IN LIVER DISEASE

Global viscoelastic (VE) tests provide a more physiological assessment of coagulation and should be considered to guide blood transfusion requirements in LT and other major surgeries. Its application in patients with ACLF or in a critical care setting requires more data. VE tests, which include TEG, ROTEM, and Sonoclot, offer a means of assessing the activity of procoagulant and anticoagulant pathways, hyperfibrinolysis, and excessive clot lysis.^{3,33} A major fallacy in the interpretation of these tests is that they are *in vitro* assays and cannot assess the *in vivo* hemostatic milieu of the endothelium, tissue factor, and portal pressure and flow. Assessment of clot formation can be performed in 10 to 20 minutes; however, assessment of clot lysis takes 30 to 60 minutes.³⁴ VE testing in LT has been described in one randomized controlled trial and several retrospective studies on the basis for which VE test guided transfusion algorithms are recommended by major guidelines.³⁵

HYPERCOAGULABILITY ASSESSMENT IN LIVER DISEASE

Patients with cirrhosis are at an increased risk (0.5%-2%) for venous thromboembolism.³⁶ Rates of PVT have been reported as approximately 8% per year, with morbidity and mortality at 1 year impacted by prophylactic anticoagulation.³⁷ Anticoagulation demonstrates the most utility in patients with more extensive portal vein and mesenteric thrombosis in the absence of other risk factors for bleeding. Low-molecular-weight heparin (LMWH) does not appear to increase risk for variceal bleeding and is likely the safest choice.³⁸ Thromboembolic events occur in 2% to 5% of patients without cirrhosis who receive thrombopoietin receptor agonists, and the risk is apparently higher in patients with cirrhosis.³¹

COMMON CLINICAL SCENARIOS IN LIVER DISEASE

Acute variceal bleeding and portal hypertensive bleeding from gastropathy, vascular ectasia, and colopathy are often challenging in patients with cirrhosis.¹¹ Another scenario encountered by the clinician is risk assessment and prophylactic correction of coagulation defects before an invasive procedure.¹⁹ Conversely, patients with cirrhosis and those with liver failure may have bleeding at one site and thrombosis at another. An example of this situation is a patient with acute PVT with extension of thrombus to the superior mesenteric vein causing mesenteric ischemia in the setting of high-grade esophageal varices. Tables 3 and 4 show these clinical scenarios and the clinical strategy to manage the condition.

DYNAMIC COAGULATION PROFILE IN SIRS AND SEPSIS

Most coagulation enzymes are also components of the inflammatory cascade. The cell-based model suggests that coagulation is one aspect of systemic inflammation. The interaction of platelets and inflammatory cytokines released from damaged endothelium triggers the adhesion of neutrophils and macrophages, and activation of inflammation. Hence patients with cirrhosis or ACLF with SIRS or sepsis are bound to have coagulation defects with an anticoagulant or procoagulant or mixed phenotype. The changes in the coagulation profile are dynamic, with resolution of global coagulation defects once sepsis resolves. Recurrence of a coagulation defect, such as platelet dysfunction, or occurrence of a new defect, such as fibrinolysis, is seen with new-onset sepsis. Prior studies have described reduced coagulation FV, FVII, FIX, FX, FXI, and prothrombin, and increased FVIII activity in cirrhosis.^{3,7-9} Two other divergent agents of the cascade affected by sepsis are tissue plasminogen activator (tPA) and plasminogen activator inhibitor (PAI). In liver failure, tPA is elevated because of endothelial activation and reduced hepatic clearance. Levels of PAI, a blocker of fibrinolysis, are also increased, but to a lesser extent than tPA.¹² Low fibrinogen is associated with a 29% increase in mortality for every 1-g/L (100-mg/dL) reduction in decompensated cirrhosis.³⁹ Hypofibrinogenemia predisposes to increased risk for bleeding after prophylactic endoscopic variceal band ligation.⁴⁰ Dysfibrinogenemia (i.e., an altered fibrinogen molecule) may cause decreased permeability of the formed clot compared with

TABLE 3. CLINICAL SCENARIOS RELATING TO COAGULATION IN CIRRHOSIS AND LIVER FAILURE

Clinical Scenario	Management	Coagulation Tests for Diagnosis and Management	Therapeutic Strategy	Pitfalls
Acute variceal bleeding	Timely endoscopic management by variceal ligation, sclerotherapy, glue injection, etc. Resuscitate as needed without increasing portal pressure.	SCT: INR PLT Fibrinogen level Global tests: TEG and ROTEM are not validated	Maintain hemoglobin target of ~7 g/dL in all patients and 8 g/dL in patients with cardiac disease Platelet transfusion target of at least $56 \times 10^9/L$ Fibrinogen level >100 mg/dL	Failure to control bleed caused by mucosal hemorrhage and early rebleeding due to high pressures Role of TIPS Role of self-expanding metal stents is of utility
Invasive procedures	Risk-to-benefit ratio and ability of the center to control bleeding should be ascertained.	INR and PLT are not a measure of coagulation risk	In high-risk procedures, PLTs of $50-60 \times 10^9/L$ are to be maintained	FFP use in the prescribed dosage of 20-40 mL/kg is unlikely to correct INR and may lead to volume overload
LT	Blood transfusion requirements are guided by the stage of the procedure with lower central venous pressures maintained during liver dissection.	Global tests: TEG and ROTEM are not validated Global tests: TEG and ROTEM are validated to guide blood component therapy	In very high-risk procedures, a PLI of $100 \times 10^9/L$ is to be targeted TEG-based correction algorithm for platelets and clotting factors in the form of FFP and cryoprecipitate	POC-based correction is validated
PVT	Anticoagulation is indicated in Yerdel grades ≥ 2 , especially in patients who are planned for LT.	Assessment of secondary fibrinolysis is essential Systemic anticoagulation with therapeutic LMWH or fondaparinux	High INR does not indicate an auto-anticoagulant state, and patients will still require anticoagulation	Patients with cavernomas are less likely to benefit
Deep vein thrombosis	In acute PVT and SMVT, catheter-mediated thrombolysis may be indicated. Consider prophylaxis in all in-hospital patients with cirrhosis as per risk scores.	Oral anticoagulants to be used include VKAs and DOACs Prophylaxis should be with LMWH or fondaparinux	Variceal eradication is necessary	Conflicting data are available on primary prophylaxis of PVT Heparin-induced thrombocytopenia and appearance of new-onset clinically evident bleeding at a remote site is challenging
Systemic thromboembolism	Anticoagulation is indicated to prevent further embolism. Assessment of variceal status and likelihood of variceal bleeding is to be weighed.	Therapeutic anticoagulation should be with oral anticoagulants INR and DIC workup is essential to guide therapy	Anticoagulation is considered for at least 6 months or lifelong in cases of multiple episodes Role of markers of fibrinolysis and DIC	
ALF	Even with increased INR, the patient tends not to bleed because of "rebalanced coagulation."	aPTT target of 2-3 times ULN measured q4-6 hours is needed Prophylactic correction of INR is NOT recommended in the absence of clinically evident bleeding	Oral anticoagulation using VKAs or DOACs may be used for secondary prophylaxis Global coagulation tests, such as TEG, often show a procoagulant state	Invasive monitoring of intracranial pressure is rarely performed, which may require correction of INR in the absence of bleeding
ACLF	Patients have coagulation failure in addition to portal hypertension and presence of renal dysfunction. If postprocedural bleeding occurs in mucosal sites or from puncture wounds, hyperfibrinolysis should be considered.	INR >1.5 is part of the definition of ACLF	INR, PLT, and fibrinogen levels are to be corrected in clinical bleeding with the same goals as cirrhosis	Role of rFVIIa (40 µg/kg) for monitor placement has been described Global tests, such as TEG and ROTEM, may have some use in managing clinical bleeding but need further validation

Abbreviations: ACLF, acute-on-chronic liver failure; ALF, acute liver failure; aPTT, activated partial thromboplastin time; DIC, disseminated intravascular coagulation; DOACs, direct oral anticoagulants; FFP, fresh frozen plasma; Hb, hemoglobin; INR, international normalized ratio; LMWH, low-molecular-weight heparin; POC, point of care; PVT, portal vein thrombosis; rF, recombinant factor; ROTEM, rotational thromboelastometry; SCT, standard coagulation test; SMVT, superior mesenteric vein thrombosis; TEG, thromboelastography; TIPS, transjugular intrahepatic portosystemic shunt; ULN, upper limit of normal.

TABLE 4. THERAPEUTIC AGENTS FOR BLEEDING OR THROMBOSIS IN LIVER DISEASE

Therapeutic Agent	Dose and Route of Administration	Mechanism of Action	Advantages	Disadvantages	Remarks
Agents for Hemostasis					
FFP	15-30 mL/kg IV	Partial correction of deficient coagulation factors	Cost-effective and available at all blood banks as a product of fractionating whole blood	Volume overload	At least 10 mL/kg is needed to effectively raise plasma protein value
Cryoprecipitate	10-20 mL plasma IV	FVIII, FXIII, fibrinogen, fibronectin, and vWF replacement	Same as FFP and provides large dose of fibrinogen in patients with DIC	Transfusion reactions	Low volume and cost-effective
Single-donor apheresis platelet unit	Platelets in a dose of 8000-12,000 cells/unit IV	Platelet boost to improve number and function	No cross-match needed	Transfusion reactions and TRALI	Use of random donor platelet units is not recommended in cirrhosis
PCCs	Typical dose is 25-30 units of FX per kg body weight IV	Three-factor (II, IX, X) and four-factor (II, IX, X, with the addition of FVII) products	Small volume	Expensive	Not widely available
rFVIIa	40-90 µg/kg IV	Also have proteins C and S	Concentrated volume	Risk for thrombosis	Use of random donor platelet units is not recommended in cirrhosis
		FVII as a driver of the coagulation cascade	Has been used in cirrhosis with good results	Expensive	Not widely available
			Concentrated volume	Risk for thrombosis: cerebrovascular and myocardial infarction and DIC	Used in acute variceal bleeding, but no mortality benefit was demonstrated
Thrombopoietin Agonists					
Eltrombopag	25-75 mg qd po × 10-20 days	Increases PLT	Expensive	Risk for PVT	Newer congeners, avatrombopag and lusutrombopag, may pose lower risk for development of PVT
			Short-lived effect		Tests of fibrinolysis are still unclear
Antifibrinolytic Agents					
Tranexamic acid	500 mg q6 hours IV or po	Disrupts the plasminogen/plasmin interaction with fibrin preventing lysis of clot	Relatively inexpensive	Doubtful risk for thrombosis or urinary clots if used in hematuria	Efficacy in cirrhosis with fibrinolysis is not yet established
EACA	Initial: 4-5 g IV/po during first hour Maintenance: 1-4 g po/IV q4-8 hours	Inhibits fibrinolysis through inhibition of plasminogen binding to fibrin and conversion to plasmin, which causes fibrinolysis	Effective topically		
Anticoagulation in Liver Disease					
Unfractionated heparin	IV or SQ injection	Inactivates FXa and thrombin	Not well studied in cirrhosis but commonly used	Needs frequent monitoring of aPTT	Reversal agent available and safe in renal disease
LMWH (enoxaparin, dalteparin)	SQ injection	Inactivates FXa and thrombin	Monitoring not needed	Extensively used in cirrhosis but expensive and possibility of heparin-induced thrombocytopenia	Reversal agent not available
VKA	po	Inhibits vitamin K-dependent FII, FVII, FIX, and FX	Cost-effective	Dosing requirement varies	Unclear dosing variables
DOACs	po	Direct thrombin inhibition by drugs such as dabigatran	INR is standard monitoring	Role in PVT	Dietary restrictions
FXa inhibitors	SQ injection (fondaparinux 7.5 mg qd) po (rivaroxaban 10 mg qd)	Direct inhibition of FXa	Oral drugs	Role in heparin-induced thrombocytopenia	CYP P450 interactions
			po or SQ administration		Monitoring and dosing protocols in cirrhosis need validation

Abbreviations: aPTT, activated partial thromboplastin time; CYP P450, cytochrome P450; DIC, disseminated intravascular coagulation; EACA, ε-aminocaproic acid; F, factor; FFP, fresh frozen plasma; INR, international normalized ratio; IV, intravenous; LMWH, low-molecular-weight heparin; PCCs, prothrombin complex concentrates; PVT, portal vein thrombosis; rF, recombinant factor; SQ, subcutaneous; TRALI, transfusion-related acute lung injury; VWF, von Willebrand factor.

controls and may even confer hypercoagulable features, causing scenarios such as portal vein or deep vein thromboses.⁴¹

The development of sepsis in patients with cirrhosis changes the coagulation factors and cellular elements, including endothelium, macrophages, and lymphocytes. Endotoxins inhibit PF by increased prostacyclin and nitric oxide production.¹² SIRS and sepsis trigger the release of endogenous heparinoids, or a heparin-like effect, caused by small endothelium/mast cell-derived glycosaminoglycans, which can be detected on heparinase-treated VE assays. One study demonstrated the presence of heparan sulphate-like molecules in patients with variceal bleeding.⁴² Senzolo et al.³³ showed that these glycosaminoglycans affect hemostasis in patients with cirrhosis and bacterial infection as seen by a hypocoagulable

native thromboelastography (TEG). This can be corrected by using heparinase-treated TEG. Bacterial sepsis predisposes patients to new bleeding events and results in difficult-to-control bleeds or endoscopic treatment failure.¹¹ Endotoxin also activates the clotting cascade and results in disseminated intravascular coagulation (DIC). After appropriate treatment, the coagulation profile resolves and endogenous heparinoids are cleared, emphasizing the complementary association of the coagulation cascade and the inflammatory pathways.¹⁷

LABORATORY TESTING IN LIVER DISEASE

SCTs, such as INR and PT, are determined by liver-synthesized coagulation FI, FII, FV, FVII, and FX. The anti-coagulation system (low protein C and protein S) is not tested, which may result in a hypercoagulable state despite


FIG 2 Coagulation correction algorithm using global coagulation tests.

increased INR in ACLF. Effective thrombin generation potential has been demonstrated to be normal in cirrhosis, AD of cirrhosis, and even uncomplicated ACLF because of higher levels of FVIII and low protein C.^{3,7-9} VE assays provide a physiological estimation of coagulation; however, their use has been validated only in the setting of LT. New research has demonstrated the use of global coagulation tests, such as TEG, ROTEM, and Sonoclot, in the setting of cirrhosis, ALF, AD of cirrhosis, and ACLF; thus, VE tests are finding mention in newer guidelines for hemostatic resuscitation, but more validation is needed before they can be recommended as standard of care.^{3,5-9,25} Table 4 shows the therapeutic targets and agents available for evidence-based management of bleeding or thromboses in liver disease. Figure 2 shows a clinical decision algorithm to determine coagulation correction targets using VE assays and SCTs.

Key elements of the whole blood TEG include the reaction time (R time), which reflects the quantity of available factors; clot formation time (K); alpha angle (α) reflecting the rate of clot formation and indirectly indicating fibrinogen levels; maximum amplitude (MA), which is an indicator of platelet activity; and finally, a measure of clot lysis.³ Thrombin generation assays provide the rate of thrombin formation in the presence of a triggering agent, such as phospholipid, a tissue factor that indicates a preserved coagulation milieu in cirrhosis or liver failure.¹³

SUMMARY AND FUTURE DIRECTIONS

The dogmatic practice of conventional testing and treatment of coagulation defects is being challenged in light of new clinical evidence of the cell-based hemostasis model in cirrhosis and liver failure. The use of global coagulation tests has provided more credible data for practical use of blood component therapy. Appreciation of the role of fibrinolysis in bleeding and effective use of antifibrinolytic agents to control mucosal bleeding are of interest in clinical practice, reducing the need for component transfusions. The pitfalls of the use of INR to guide transfusion strategies are well recognized. Hypercoagulable states in the setting of cirrhosis and liver failure are often undiagnosed. Prospective studies are needed to examine the utility of current risk models in patients with liver disease with a goal of improving knowledge of the incidence and natural history of venous thromboembolism.

CORRESPONDENCE

Shiv Kumar Sarin, M.D., D.M., Department of Hepatology, Institute of Liver and Biliary Sciences, D1 Vasant Kunj, New Delhi 110070, India.
E-mail: sksarin@ilbs.in

REFERENCES

- 1) Lisman T, Caldwell SH, Burroughs AK, et al.; Coagulation in Liver Disease Study Group. Hemostasis and thrombosis in patients with liver disease: the ups and downs. *J Hepatol* 2010;53:362-371.
- 2) Tripodi A, Primignani M, Mannucci PM, et al. Changing concepts of cirrhotic coagulopathy. *Am J Gastroenterol* 2017;112:274-281.
- 3) Premkumar M, Saxena P, Rangegowda D, et al. Coagulation failure is associated with bleeding events and clinical outcome during systemic inflammatory response and sepsis in acute-on-chronic liver failure: an observational cohort study. *Liver Int* 2019;39:694-704.
- 4) Northup PG, Caldwell SH. Coagulation in liver disease: a guide for the clinician. *Clin Gastroenterol Hepatol* 2013;11:1064-1074.
- 5) Agarwal B, Wright G, Gatt A, et al. Evaluation of coagulation abnormalities in acute liver failure. *J Hepatol* 2012;57:780-786.
- 6) Sarin SK, Choudhury A, Sharma MK, et al.; APASL ACLF Research Consortium (AARC) for APASL ACLF working Party. Acute-on-chronic liver failure: consensus recommendations of the Asian Pacific association for the study of the liver (APASL): an update [corrected]. *Hepatol Int* 2019;13:826-828. *Hepatol Int* 2019;13:353-390.
- 7) Blasi A, Calvo A, Prado V, et al. Coagulation failure in patients with acute-on-chronic liver failure (ACLF) and decompensated cirrhosis: beyond INR. *Hepatology* 2018;68:2325-2337.
- 8) Fisher C, Patel VC, Stoy SH, et al. Balanced haemostasis with both hypo- and hyper-coagulable features in critically ill patients with acute-on-chronic-liver failure. *J Crit Care*. 2018;43:54-60. doi:10.1016/j.jcrc.2017.07.053
- 9) Bedreli S, Sowa JP, Gerken G, et al. Management of acute-on-chronic liver failure: rotational thromboelastometry may reduce substitution of coagulation factors in liver cirrhosis. *Gut* 2016;65:357-358.
- 10) Case JJ, Khan N, Delrahim M, et al. Association of massive transfusion for resuscitation in gastrointestinal bleeding with transfusion-related acute lung injury. *Indian J Crit Care Med* 2017;21:506-513.
- 11) Kumar M, Ahmad J, Maiwall R, et al. Thromboelastography-guided blood component use in patients with cirrhosis with nonvariceal bleeding: a randomized controlled trial. *Hepatology* 2020;71:235-246.
- 12) Tripodi A. Hemostasis abnormalities in cirrhosis. *Curr Opin Hematol* 2015;22:406-412.
- 13) Tripodi A, Primignani M, Chantarangkul V, et al. Thrombin generation in patients with cirrhosis: the role of platelets. *Hepatology* 2006;44:440-445.
- 14) Tripodi A, Salerno F, Chantarangkul V, et al. Evidence of normal thrombin generation in cirrhosis despite abnormal conventional coagulation tests. *Hepatology* 2005;41:553-558.

- 15) Silva MA, Jambulingam PS, Gunson BK, et al. Hepatic artery thrombosis following orthotopic liver transplantation: a 10-year experience from a single centre in the United Kingdom. *Liver Transpl* 2006;12:146-151.
- 16) Salami A, Qureshi W, Kuriakose P, et al. Frequency and predictors of venous thromboembolism in orthotopic liver transplant recipients: a single-center retrospective review. *Transplant Proc* 2013;45:315-319.
- 17) Premkumar M, Saxena P, Mirza P, et al. Heparin like effect increases the risk of mortality and bleeding in patients with severe alcoholic hepatitis. *Clin Gastroenterol Hepatol* 2020;18:486-495.
- 18) Fayed N, Mourad W, Yassen K, et al. Preoperative thromboelastometry as a predictor of transfusion requirements during adult living donor liver transplantation. *Transfus Med Hemother* 2015;42:99-108.
- 19) De Pietri L, Bianchini M, Montalti R, et al. Thrombelastography-guided blood product use before invasive procedures in cirrhosis with severe coagulopathy: a randomized, controlled trial. *Hepatology* 2016;63:566-573.
- 20) Stravitz RT, Lisman T, Luketic VA, et al. Minimal effects of acute liver injury/acute liver failure on hemostasis as assessed by thromboelastography. *J Hepatol* 2012;56:129-136.
- 21) Goyal S, Jadaun S, Kedia S, et al. Thromboelastography parameters in patients with acute on chronic liver failure. *Ann Hepatol* 2018;17:1042-1051.
- 22) Papatheodoridis GV, Patch D, Webster GJ, et al. Infection and hemostasis in decompensated cirrhosis: a prospective study using thrombelastography. *Hepatology* 1999;29:1085-1090.
- 23) Lloyd-Donald P, Vasudevan A, Angus P, et al. Coagulation in acutely ill patients with severe chronic liver disease: insights from thromboelastography. *J Crit Care* 2017;38:215-224.
- 24) Fayed N, Mourad W, Yassen K, et al. Preoperative thromboelastometry as a predictor of transfusion requirements during adult living donor liver transplantation. *Transfus Med Hemother* 2015;42:99-108.
- 25) Intagliata NM, Argo CK, Stine JG, et al.; Faculty of the 7th International Coagulation in Liver Disease. Concepts and controversies in haemostasis and thrombosis associated with liver disease: Proceedings of the 7th International Coagulation in Liver Disease Conference. *Thromb Haemost* 2018;118:1491-1506.
- 26) Kozek-Langenecker SA, Afshari A, Albaladejo P, et al. Management of severe perioperative bleeding: guidelines from the European Society of Anaesthesiology. *Eur J Anaesthesiol* 2013;30:270-382.
- 27) Stravitz RT. Algorithms for managing coagulation disorders in liver disease. *Hepatol Int* 2018;12:390-401.
- 28) Zimmon DS, Kessler RE. The portal pressure-blood volume relationship in cirrhosis. *Gut* 1974;15:99-101.
- 29) Kwon JO, MacLaren R. Comparison of fresh-frozen plasma, four-factor prothrombin complex concentrates, and recombinant factor VIIa to facilitate procedures in critically ill patients with coagulopathy from liver disease: A retrospective cohort study. *Pharmacotherapy* 2016;36:1047-1054.
- 30) Afdhal NH, Giannini EG, Tayyab G, et al.; ELEVATE Study Group. Eltrombopag before procedures in patients with cirrhosis and thrombocytopenia. *N Engl J Med* 2012;367:716-724.
- 31) Michelson AD, Smolensky Koganov E, Forde EE, et al. Avatrombopag increases platelet count but not platelet activation in patients with thrombocytopenia resulting from liver disease. *J Thromb Haemost* 2018;16:2515-2519.
- 32) Terrault NA, Hassanein T, Howell CD, et al. Phase II study of avatrombopag in thrombocytopenic patients with cirrhosis undergoing an elective procedure. *J Hepatol* 2014;61:1253-1259.
- 33) Senzolo M, Coppell J, Cholangitas E, et al. The effects of glycosaminoglycans on coagulation: a thromboelastographic study. *Blood Coagul Fibrinolysis* 2007;18:227-236.
- 34) Krzanicki D, Sugavanam A, Mallett S. Intraoperative hypercoagulability during liver transplantation as demonstrated by thromboelastography. *Liver Transpl* 2013;19:852-861.
- 35) Northup PG, Friedman LS, Kamath PS. AGA clinical practice update on surgical risk assessment and perioperative management in cirrhosis: expert review. *Clin Gastroenterol Hepatol* 2019;17:595-606.
- 36) Gulley D, Teal E, Suvannasankha A, et al. Deep vein thrombosis and pulmonary embolism in cirrhosis patients. *Dig Dis Sci* 2008;53:3012-3017.
- 37) Faccia M, Ainora ME, Ponziani FR, et al. Portal vein thrombosis in cirrhosis: why a well-known complication is still matter of debate. *World J Gastroenterol* 2019;25:4437-4451.
- 38) Villa E, Cammà C, Marietta M, et al. Enoxaparin prevents portal vein thrombosis and liver decompensation in patients with advanced cirrhosis. *Gastroenterology* 2012;143:1253-1260.e1-e4.
- 39) Desborough MJ, Kahan BC, Stanworth SJ, et al. Fibrinogen as an independent predictor of mortality in decompensated cirrhosis and bleeding. *Hepatology* 2017;65:1079-1080.
- 40) Giannini EG, Giambruno E, Brunacci M, et al. Low fibrinogen levels are associated with bleeding after varices ligation in thrombocytopenic cirrhotic patients. *Ann Hepatol* 2018;17:830-835.
- 41) Hugenholtz GC, Macrae F, Adelmeijer J, et al. Procoagulant changes in fibrin clot structure in patients with cirrhosis are associated with oxidative modifications of fibrinogen. *J Thromb Haemost* 2016;14:1054-1066.
- 42) Montalto P, Vlachogiannakos J, Cox DJ, et al. Bacterial infection in cirrhosis impairs coagulation by a heparin effect: a prospective study. *J Hepatol* 2002;37:463-447.

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IMPROVE GUANG ZHOU IMPROVE MEDICAL INSTRUMENTS CO.,LTD