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IMPROVE REVIEW

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Thromboelastography

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- SHAYDAKOV M E, SIGMON D F, BLEBEA J. Thromboelastography [M]. StatPearls. Treasure Island (FL) ineligible companies. Disclosure: David Sigmon declares no relevant financial relationships with ineligible companies. Disclosure: John Blebea declares no relevant financial relationships with ineligible companies. 2023.
- HARTERT H. [Blood clotting studies with Thrombus stressography; a new Investigation procedure] [J]. Klin Wochenschr, 1948, 26(37-38): 577-83.
- WALSH M, FRITZ S, HAKE D, et al. Targeted Thromboelastographic (TEG) Blood Component and Pharmacologic Hemostatic Therapy in Traumatic and Acquired Coagulopathy [J]. Curr Drug Targets, 2016, 17(8): 954-70.
- KANG Y G, MARTIN D J, MARQUEZ J, et al. Intraoperative changes in blood coagulation and thromboelastographic monitoring in liver transplantation [J]. Anesthesia and Analgesia, 1985, 64(9): 888-96.
- SHORE-LESSERSON L, MANSPEIZER H E, DEPERIO M, et al. Thromboelastography-guided transfusion algorithm reduces transfusions in complex cardiac surgery [J]. Anesthesia and Analgesia, 1999, 88(2): 312-9.
- 赵晖, 李斌, 张元珍, et al. 血栓弹力图在实验分析前影响因素的探讨 [J]. 中国卫生检验杂志, 2015, 25(13): 2139-40.
- 王海, 邱平, 齐晓伟. 血栓弹力图实验操作影响因素的探讨 [J]. 血栓与止血学, 2011, 17(06): 272-4.
- LIPPI G, SALVAGNO G L, MONTAGNANA M, et al. Quality standards for sample collection in coagulation testing [J]. Semin Thromb Hemost, 2012, 38(6): 565-75.
- WALENGA J M, JACKSON C M, KESSLER C M. Low molecular weight heparins differ substantially: impact on developing biosimilar drugs [J]. Seminars In Thrombosis and Hemostasis, 2011, 37(3): 322-7.
- TUFANO A, COPPOLA A, CERBONE A M, et al. Preventing postsurgical venous thromboembolism: pharmacological approaches [J]. Seminars In Thrombosis and Hemostasis, 2011, 37(3): 252-66.
- FAVALORO E J, LIPPI G, KOUTTS J. Laboratory testing of anticoagulants: the present and the future [J]. Pathology, 2011, 43(7): 682-92.
- FAVALORO E J, LIPPI G. Laboratory testing and/or monitoring of the new oral anticoagulants/antithrombotics: for and against? [J]. Clin Chem Lab Med, 2011, 49(5): 755-7.
- LIPPI G, FAVALORO E J, SALVAGNO G L, et al. Laboratory assessment and perioperative management of patients on antiplatelet therapy: from the bench to the bedside [J]. Clin Chim Acta, 2009, 405(1-2):
- FAVALORO E J, LIPPI G. Laboratory reporting of hemostasis assays: the final post-analytical opportunity to reduce errors of clinical diagnosis in hemostasis? [J]. Clin Chem Lab Med, 2010, 48(3): 309-21.
- LIPPI G, SALVAGNO G L, MONTAGNANA M, et al. Influence of the needle bore size used for collecting venous blood samples on routine clinical chemistry testing [J]. Clin Chem Lab Med, 2006, 44(8): 1009-14.
- LIPPI G, SALVAGNO G L, MONTAGNANA M, et al. Influence of the needle bore size on platelet count and routine coagulation testing [J]. Blood Coagul Fibrinolysis, 2006, 17(7): 557-61.
- GOSSELIN R C, JANATPOUR K, LARKIN E C, et al. Comparison of samples obtained from 3.2% sodium citrate glass and two 3.2% sodium citrate plastic blood collection tubes used in coagulation testing [J]. Am J Clin Pathol, 2004, 122(6): 843-8.
- FIEBIG E W, ETZELL J E, NG V L. Clinically relevant differences in prothrombin time and INR values related to blood sample collection in plastic vs glass tubes [J]. Am J Clin Pathol, 2005, 124(6): 902-9.
- KRATZ A, STANGANELLI N, VAN COTT E M. A comparison of glass and plastic blood collection tubes for routine and specialized coagulation assays: a comprehensive study [J]. Arch Pathol Lab Med, 2006, 130(1): 39-44.
- TOULON P, AJZENBERG N, SMAHI M, et al. A new plastic collection tube made of polyethylene terephthalate is suitable for monitoring traditional anticoagulant therapy (oral anticoagulant, unfractionated heparin, and low molecular weight heparin) [J]. Thromb Res, 2007, 119(2): 135-43.
- 李贤见, 王长本, 唐李. 不同品牌采血管和标本放置时间对血栓弹力图检测结果的影响分析 [J]. 标记免疫分析与临床, 2019, 26(10): 1770-4.
- PRESTON F E, LIPPI G, FAVALORO E J, et al. Quality issues in laboratory haemostasis [J]. Haemophilia, 2010, 16 Suppl 5(93-9).
- LIPPI G, SALVAGNO G L, MONTAGNANA M, et al. Quality standards for sample collection in coagulation testing [J]. Seminars In Thrombosis and Hemostasis, 2012, 38(6): 565-75.
- ADCOCK FUNK D M, LIPPI G, FAVALORO E J. Quality standards for sample processing, transportation, and storage in hemostasis testing [J]. Semin Thromb Hemost, 2012, 38(6): 576-85.
- HEIL W, GRUNEWALD R, AMEND M, et al. Influence of time and temperature on coagulation analytes in stored plasma [J]. Clin Chem Lab Med, 1998, 36(7): 459-62.
- JÖNSSON RYLANDER A C, LINDGREN A, DEINUM J, et al. Fibrinolysis inhibitors in plaque stability: a morphological association of PAI-1 and TAFI in advanced carotid plaque [J]. J Thromb Haemost, 2017, 15(4): 758-69.
- KATTLOVE H E, ALEXANDER B. The effect of cold on platelets. I. Cold-induced platelet aggregation [J]. Blood, 1971, 38(1): 39-48.
- KRISTOFFERSEN A H, STAVELIN A V, AJZNER E, et al. Pre-analytical practices for routine coagulation tests in European laboratories. A collaborative study from the European Organisation for External Quality Assurance Providers in Laboratory Medicine (EQALM) [J]. Clin Chem Lab Med, 2019, 57(10): 1511-21.
- 张海霞. 血标本放置时间对血液分析仪血小板计数的影响 [J]. 上海医学检验杂志, 2002, 02(02): 109-10.

30. LEMON A V, GODDARD A, HOOLJBERG E H. Effects of storage time and temperature on thromboelastographic analysis in dogs and horses [J]. *Vet Clin Pathol*, 2021, 50(1):
31. VIG S, CHITOLIE A, BEVAN D H, et al. Thromboelastography: a reliable test? [J]. *Blood Coagul Fibrinolysis*, 2001, 12(7): 555-61.
32. QUARTERMAN C, SHAW M, JOHNSON I, et al. Intra- and inter-centre standardisation of thromboelastography (TEG®) [J]. *Anaesthesia*, 2014, 69(8): 883-90.
33. 陈勇, 凌建, 徐容. 影响凝血四项检测结果异常原因的分析 [J]. *中国实验诊断学*, 2017, 21(09): 1523-7.
34. 马跃飞, 林锦骝, 吴文楠. 血栓弹力图分析前影响因素的探讨 [J]. *中国卫生检验杂志*, 2020.20(09): 1127-30.
35. CLINICAL, LABORATORY STANDARDS I. Tubes and Additives for Venous Blood Specimen Collection; Approved Standard. 5th ed. CLSI Document H1-A5 [M]. Wayne, PA: Clinical Laboratory Standards Institute, 2003.
36. CLINICAL, LABORATORY STANDARDS I. Platelet Function Testing by Aggregometry; Approved Guideline. CLSI Document H58-A [M]. Wayne, PA: Clinical Laboratory Standards Institute, 2008.
37. CATTANEO M. Light transmission aggregometry and ATP release for the diagnostic assessment of platelet function [J]. *Seminars In Thrombosis and Hemostasis*, 2009, 35(2): 158-67.
38. DYSZKIEWICZ-KORPANTY A, QUINTON R, JASSINE J, et al. The effect of a pneumatic tube transport system on PFA-100 trade mark closure time and whole blood platelet aggregation [J]. *J Thromb Haemost*, 2004, 2(354-6).
39. PLEBANI M, ZANINOTTO M. Pneumatic tube delivery systems for patient samples: evidence of quality and quality of evidence [J]. *Clin Chem Lab Med*, 2011, 49(8): 1245-6.
40. KRATZ A, SALEM R O, VAN COTT E M. Effects of a pneumatic tube system on routine and novel hematology and coagulation parameters in healthy volunteers [J]. *Arch Pathol Lab Med*, 2007, 131(2): 293-6.
41. WALLIN O, SÖDERBERG J, GRANKVIST K, et al. Preanalytical effects of pneumatic tube transport on routine haematology, coagulation parameters, platelet function and global coagulation [J]. *Clin Chem Lab Med*, 2008, 46(10): 1443-9.
42. MONTARULI B, GUIOTTO C, COSSEDDU D. Influence of hemolysis, icterus and lipemia on coagulation tests as performed on Cobas t511 new analyzer [J]. *Blood Coagul Fibrinolysis*, 2020, 31(1): 48-54.
43. FAVALORO E J, LIPPI G. Recent Advances in Mainstream Hemostasis Diagnostics and Coagulation Testing [J]. *Seminars In Thrombosis and Hemostasis*, 2019, 45(3): 228-46.
44. LIPPI G, BETSOU F, CADAMURO J, et al. Preanalytical challenges - time for solutions [J]. *Clin Chem Lab Med*, 2019, 57(7): 974-81.
45. WHITTEN C W, GREILICH P E. Thromboelastography: past, present, and future [J]. *Anesthesiology*, 2000, 92(5): 1223-5.
46. TANG N, JIN X, SUN Z, et al. Effects of hemolysis and lipemia interference on kaolin-activated thromboelastography, and comparison with conventional coagulation tests [J]. *Scand J Clin Lab. Invest*, 2017, 77(2):
47. LIPPI G, BLANCKAERT N, BONINI P. Haemolysis: an overview of the leading cause of unsuitable specimens in clinical laboratories [J]. *Clin Chem Lab Med*, 2008, 46(6): 764-72.
48. LIPPI G, PLEBANI M, DI SOMMA S, et al. Hemolyzed specimens: a major challenge for emergency departments and clinical laboratories [J]. *Crit Rev Clin Lab Sci*, 2011, 48(3): 143-53.
49. LIPPI G, MUSA R, AVANZINI P, et al. Influence of in vitro hemolysis on hematological testing on Advia 2120 [J]. *Int J Lab Hematol*, 2012, 34(2): 179-84.
50. CLINICAL LABORATORY STANDARDS I. Platelet Function Testing by Aggregometry; Approved Guideline. CLSI document H58-A [M]. Wayne, PA: Clinical and Laboratory Standards Institute, 2008.
51. KITCHEN S, ADCOCK D M, DAUER R, et al. International Council for Standardization in Haematology (ICSH) recommendations for processing of blood samples for coagulation testing [J]. *Int J Lab Hematol*, 2021, 43(6): 1272-83.
52. CLINICAL, LABORATORY STANDARDS I. Quality Management System: A Model for Laboratory Services; Approved Guideline. 4th ed. CLSI Document GP26-A4 [M]. Wayne, PA: Clinical and Laboratory Standards Institute, 2011.
53. CLINICAL, LABORATORY STANDARDS I. Training and Competence Assessment; Approved Guideline. 3rd ed. CLSI Document GP21-A3 [M]. Wayne, PA: Clinical and Laboratory Standards Institute, 2009.
54. HAMMERLING J A. A review of medical errors in laboratory diagnostics and where we are today [J]. *Lab Med*, 2012, 43(41-4).
55. FAVALORO E J, LIPPI G C, ADCOCK D M. Preanalytical and postanalytical variables: the leading causes of diagnostic error in hemostasis? [J]. *Semin Thromb Hemost*, 2008, 34(7): 612-34.
56. FAVALORO E J, ADCOCK D M, LIPPI G C. Preanalytical variables in coagulation testing associated with diagnostic errors in hemostasis [J]. *Lab Med*, 2012, 43(2): 1-10.
57. LIPPI G, CHANCE J J, CHURCH S. Preanalytical quality improvement: from dream to reality [J]. *Clin Chem Lab Med*, 2011, 49(7): 1113-26.
58. SELBY R, DUFFY E R, BERNY-LANG M A, et al. Impact of specific preclinical variables on coagulation biomarkers in cancer-associated thrombosis [J]. *Thromb Res*, 2020, 191 Suppl 1(S26-S30).
59. FAVALORO E J, LIPPI G, ADCOCK D M. Preanalytical and postanalytical variables: the leading causes of diagnostic error in hemostasis? [J]. *Seminars In Thrombosis and Hemostasis*, 2008, 34(7): 612-34
60. ZÜRCHER M, SULZER I, BARIZZI G, et al. Stability of coagulation assays performed in plasma from citrated whole blood transported at ambient temperature [J]. *Thrombosis and Haemostasis*, 2008, 99(2): 416-26.

61. LINSKENS E A, DEVREESE K M J. Pre-analytical stability of coagulation parameters in plasma stored at room temperature [J]. *Int J Lab Hematol*, 2018, 40(3): 292-303.
62. BAKER P, PLATTON S, GIBSON C, et al. Guidelines on the laboratory aspects of assays used in haemostasis and thrombosis [J]. *Br J Haematol*, 2020, 191(3): 347-62.
63. LIPPI G, PLEBANI M, FAVALORO E J. Interference in coagulation testing: focus on spurious hemolysis, icterus, and lipemia [J]. *Seminars In Thrombosis and Hemostasis*, 2013, 39(3): 258-66.
64. MAGNETTE A, CHATELAIN M, CHATELAIN B, et al. Pre-analytical issues in the haemostasis laboratory: guidance for the clinical laboratories [J]. *Thromb J*, 2016, 14(49).
65. ADCOCK D M, KRESSIN D C, MARLAR R A. Minimum specimen volume requirements for routine coagulation testing: dependence on citrate concentration [J]. *Am J Clin Pathol*, 1998, 109(5): 595-9.
66. LAGA A C, CHEVES T A, SWEENEY J D. The effect of specimen hemolysis on coagulation test results [J]. *Am J Clin Pathol*, 2006, 126(5): 748-55.
67. 赵欣宇, 宋鉴清, 张丽霞. 标本放置时间对凝血酶原时间和活化的部分凝血活酶时间检测结果的影响 [J]. *中国医科大学学报*, 2012, 41(7): 3.

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1. Lippi, Giuseppe et al. "Quality and reliability of routine coagulation testing: can we trust that sample?." *Blood coagulation & fibrinolysis : an international journal in haemostasis and thrombosis* vol. 17,7 (2006): 513-9. doi:10.1097/01.mbc.0000245290.57021.46.

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2. Tang, Ning et al. "Effects of hemolysis and lipemia interference on kaolin-activated thromboelastography, and comparison with conventional coagulation tests." *Scandinavian journal of clinical and laboratory investigation* vol. 77,2 (2017): 98-103. doi:10.1080/00365513.2016.1271906.

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
3. Cheng, Daye et al. "Establishment of thromboelastography reference intervals by indirect method and relevant factor analyses." *Journal of clinical laboratory analysis* vol. 34,6 (2020): e23224. doi:10.1002/jcla.23224.

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4. Gilman, Elizabeth A et al. "Fresh and citrated whole-blood specimens can produce different thromboelastography results in patients on extracorporeal membrane oxygenation." *American journal of clinical pathology* vol. 140,2 (2013): 165-9. doi:10.1309/AJCPYIQ9JNNSEN4Q

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5. Dias, Joao D et al. "New-Generation Thromboelastography: Comprehensive Evaluation of Citrated and Heparinized Blood Sample Storage Effect on Clot-Forming Variables." *Archives of pathology & laboratory medicine* vol. 141,4 (2017): 569-577. doi:10.5858/arpa.2016-0088-OA.

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
7. Kristoffersen, Ann Helen et al. "Pre-analytical practices for routine coagulation tests in European laboratories. A collaborative study from the European Organisation for External Quality Assurance Providers in Laboratory Medicine (EQALM)." *Clinical chemistry and laboratory medicine* vol. 57,10 (2019): 1511-1521. doi:10.1515/cclm-2019-0214.

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8. Keppel, Martin H et al. "Heparin and citrate additive carryover during blood collection." *Clinical chemistry and laboratory medicine* vol. 57,12 (2019): 1888-1896. doi:10.1515/cclm-2019-0433.

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9. Luo, Cuizhu et al. "Thromboelastography Parameters as Predictors for Long-Term Survival in Critically Ill Patients." *Clinical and applied thrombosis/hemostasis : official journal of the International Academy of Clinical and Applied Thrombosis/Hemostasis* vol. 25 (2019): 1076029619876028. doi:10.1177/1076029619876028.

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10. Kitchen, Steve et al. "International Council for Standardization in Haematology (ICSH) recommendations for processing of blood samples for coagulation testing." *International journal of laboratory hematology* vol. 43,6 (2021): 1272-1283. doi:10.1111/ijlh.13702

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12. White, H et al. "Can Thromboelastography performed on kaolin-activated citrated samples from critically ill patients provide stable and consistent parameters?." *International journal of laboratory hematology* vol. 32,2 (2010): 167-73. doi:10.1111/j.1751-553X.2009.01152.x.

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International Council for Standardization in Haematology (ICSH) recommendations for processing of blood samples for coagulation testing

Steve Kitchen¹  | Dorothy M. Adcock²  | Ray Dauer³ | Ann-Helen Kristoffersen^{4,5} | Giuseppe Lippi⁶  | Ian Mackie⁷  | Richard A. Marlar⁸  | Sukesh Nair⁹ 

¹Sheffield Haemophilia and Thrombosis Centre, Sheffield, UK

²Laboratory Corporation of America Holdings, Burlington, North Carolina, USA

³Peter MacCallum Cancer Centre, Melbourne, Victoria, Australia

⁴Department of Medical Biochemistry and Pharmacology, Haukeland University Hospital, Bergen, Norway

⁵Norwegian Organization for Quality Improvement of Laboratory Examinations (Noklus), Haraldsplass Deaconess Hospital, Bergen, Norway

⁶Section of Clinical Biochemistry, University of Verona, Verona, Italy

⁷Research Department of Haematology, University College London, London, UK

⁸Department of Pathology, University of New Mexico, Albuquerque, New Mexico, USA

⁹Christian Medical College, Vellore, India

Correspondence

Steve Kitchen, Haematology Scientist, Peter MacCallum Cancer Centre, 305 Grattan St., Melbourne, Vic., 3000.
Email: Steve.kitchen1@nhs.net

Funding information

The drafting group are grateful to the Western Norway Regional Health Authority for supporting Ann Helen Kristoffersen with postdoctoral fellowships to facilitate participation in drafting this guidance.

Abstract

This guidance document has been prepared on behalf of the International Council for Standardization in Haematology (ICSH). The aim of the document is to provide guidance and recommendations for the processing of citrated blood samples for coagulation tests in clinical laboratories in all regions of the world. The following areas are included in this document: Sample transport including use of pneumatic tubes systems; clots in citrated samples; centrifugation; primary tube storage and stability; interfering substances including haemolysis, icterus and lipaemia; secondary aliquots—transport, storage and processing; preanalytical variables for platelet function testing. The following areas are excluded from this document, but are included in an associated ICSH document addressing collection of samples for coagulation tests in clinical laboratories; ordering tests; sample collection tube and anticoagulant; preparation of the patient; sample collection device; venous stasis before sample collection; order of draw when different sample types are collected; sample labelling; blood-to-anticoagulant ratio (tube filling); influence of haematocrit. The recommendations are based on published data in peer-reviewed literature and expert opinion.

KEYWORDS

coagulation, ICSH, sample processing

1 | INTRODUCTION

Preanalytical variables affect the results of many laboratory tests including those related to haemostasis and thrombosis.^{1,2} The preanalytical phase is often the main source of inaccurate laboratory test results,³ with between one third and three quarters of laboratory errors being attributable to this phase.⁴ Many preanalytical errors are a consequence of inappropriate or problematic blood sample

processing, and there are significant cost implications of preanalytical errors as indicated by a study which estimated the average cost of a preanalytical error to be around \$200 (in 2013) in both European and North American institutions, representing an annual cost of \$1.2 million for a 650 bed hospital in the USA.⁵

Preanalytical problems can be surprisingly common occurring in around 5% of samples in one study.⁶ The analysis of unsuitable or inappropriate samples and the release of results unsafe for informing

clinical management decisions constitute important risks for patient safety. Preanalytical problems associated with sample collection are described in a separate ICSH document.⁷ The following sections describe a number of circumstances when preanalytical problems may lead to bias in test results to the extent that patient safety could be compromised.

2 | SAMPLE TRANSPORT

Fresh samples should be transported to and within the laboratory promptly at ambient (18-24°C) temperature.⁸ Storage of whole blood samples on ice for up to 4 hours affects some coagulation test results in normal subjects but not others.⁹ Transporting samples on ice or storing whole blood in a refrigerator prior to processing is not recommended for most coagulation tests due to the potential for cold precipitation of clotting factors including fibrinogen, von Willebrand (VWF) factor and factor VIII,¹⁰ as well as cold activation of factor VII and FXII with shortening of activated partial thromboplastin time (APTT) and especially prothrombin time (PT).¹¹ PT assessed with thromboplastins containing bovine tissue factor is particularly affected by cold activation of FVII.¹⁰ Laboratories must ensure that users of their services are aware of the importance of avoiding storage of whole blood samples for coagulation tests at 2-8°C. Citrated samples should not be stored at 2-8°C prior to analysis as this can lead to time-dependent loss of FVIII and VWF in whole blood^{12,13} to such an extent that normal subjects could be misclassified as having von Willebrand disease (VWD). In one study, 50% of samples stored for 3.5 hours at 2-8°C before centrifugation and analysis were misdiagnosed as VWD. Although this marked reduction could be reversed if the whole blood sample is pre-warmed prior to centrifugation,¹⁴ these problems are best avoided by transport and storage of samples at 18-24°C.

Looking to possible new transport systems in the future, preliminary data on transport of citrated blood samples by unmanned aerial drones for 6-38 minutes flights indicated no clinically relevant bias for PT or APTT.¹⁵

Issues related to transport of secondary aliquots are addressed in Section 8 below.

Transport of both primary tubes and secondary aliquots must take account of any applicable guidance or requirements from regulatory agencies on transporting biohazardous materials.

Recommendation 2.1: Whole blood samples for coagulation tests should be maintained at ambient temperature of 18-24°C during transport and storage prior to processing.

3 | PNEUMATIC TUBE SYSTEMS

Transport of coagulation samples through certain pneumatic tube systems (PTS) is acceptable for some coagulation tests depending on the conditions,¹⁶ although vibration and movement of the receptacle used for transport can lead to platelet activation in citrated blood and should not be used prior to tests of platelet function.¹⁷ PTS can

induce in vitro haemolysis as specimens are transported from inpatient wards and local outpatient departments (ODP) to the laboratory for testing.^{18,19}

While cost-effective and time efficient,²⁰ PTS can have a clinically significant effect on quality of samples, as they are exposed to different temperatures, vibrations and sudden changes in acceleration (g-force).²¹ Such physical forces contribute to the rupture of red blood cells (PTS-induced cell lysis)²² leading to an increase of intracellular constituents, including haemoglobin, in the surrounding plasma.^{23,24} Haemolysis in vitro also activates platelets with implications for coagulation laboratory tests.²⁵ Studies have indicated that acceleration forces exerted upon samples during PTS transport should be measured, to confirm a transport system's reliability and effect on sample integrity.²⁶ The relevant International Standards Organization document, ISO 15189 requires preanalytical specimen transport to be monitored,²⁷ and miniaturized data loggers have been used to measure acceleration during sample transport through PTS. Some measures of sample movement during PTS transport correlate with the rates of sample haemolysis, where it is likely that numerous changes in acceleration across multiple axes create sufficient shear stresses to induce haemolysis.^{26,28} Appropriate filling of vacant airspace in PTS carriers with packaging significantly reduces vibration and consequent haemolysis.²⁹ Each laboratory should evaluate their own PTS for induction of haemolysis in coagulation samples. Laboratories should do more than simply identify the presence of in vitro haemolysis. Using data loggers allows modifications of transport methods to be assessed, quality improvements to be made, rejection rates to be reduced and ultimately, patient care to be improved. For these reasons, the use of data loggers is recommended to assess potential impact of specimen transport via PTS.

Recommendation 3.1: Pneumatic tube transport systems should not be used to transport citrated blood samples prior to tests of platelet function.

Recommendation 3.2: Each laboratory should assess the impact of transport of citrated blood samples through their pneumatic tube system for clotting screening tests prior to any routine use of such systems.

4 | CLOTS IN CITRATED BLOOD SAMPLES

The coagulation process can be activated during collection and processing of citrated blood samples. This can lead to generation of activated clotting factors, which may result in shortening of some clotting times. The clotting process within the collection tube can lead to generation of clots which may or may not be visible on scrutiny. There can be clots inside the sedimented red cell mass after centrifugation, there can be small visible clots containing entrapped red blood cells suspended in the plasma of centrifuged samples and there can be clotting in the plasma over the top of red cell mass after centrifugation has been completed. Activated samples may contain fibrin clots which are pale in colour. Fibrin clots may be apparent after samples have gone through a freeze-thaw cycle. Visible

inspection on gentle inversion can sometimes identify the presence of clots, but it is not always possible to detect all clots. This means that assessment of results should include consideration of possible effects of activation or clotting in the sample. Formation of clots leads to consumption of factor II, factor V, factor VIII and fibrinogen from plasma. For this reason, the PT and APTT on clotted samples are usually unreadable, so that automated systems fail to record a clotting time. Fibrinogen is undetectable in fully clotted citrated samples. If there is any evidence of clotting of citrated blood samples for coagulation testing, the samples should be rejected.

Partial clotting or activation may be more difficult to identify, but may be associated with shortened APTT values. APTT may even be shorter than those observed in conditions associated with elevated clotting factors, such as third trimester pregnancy,³⁰ stress,³¹ acute phase reaction or presence of emicizumab.³² Most APTT results which are 1-3 seconds (s) below the lower limit of the reference range are likely to be genuine results caused by elevated intrinsic clotting factors, particularly FVIII. The presence of unexpected APTT results that are markedly reduced, typically for example ≥ 4 seconds below the lower limit of the reference range, should be followed by further sample scrutiny to look for clots, with consideration given to requesting a replacement sample even if clots are not apparent. The combination of very short APTT and low fibrinogen can occur in conditions associated with pathological or therapeutic fibrinolysis or fibrinogen lysis, but is more commonly an artefact of sample activation/clotting.

Recommendation 4.1: Any citrated blood samples, for coagulation tests, should be rejected if clots of any size are detected.

Recommendation 4.2: Citrated samples with APTT results 4 or more seconds below the lower limit of reference range should be scrutinized for possible presence of clots.

5 | CENTRIFUGATION

Centrifugation of citrate samples for most coagulation tests should be sufficient to deliver platelet-poor plasma with a residual platelet count of $<10 \times 10^9/L$ ⁸ which can be achieved by centrifugation at minimum of 1500 g for 15 minutes³³ or 1700 g for 10 minutes. Platelet counts of platelet-poor plasma after centrifugation should be checked before introducing new centrifuges into routine use, after any relevant centrifuge repair or recalibration, and taking account of any relevant regulatory requirements. There are some data which indicate a lack of impact of residual platelet counts up to $200 \times 10^9/L$ on PT, APTT and thrombin time when testing fresh samples.³⁴ Centrifugation at higher speeds for shorter periods at higher g-force may also be acceptable. One study validated the use of centrifugation at 11 000 g for 2 minutes for PT, APTT and fibrinogen testing,³⁵ although others reported that platelets can be released from the tube wall back into plasma where fixed angle rotors are used, when testing was delayed more than a few minutes after centrifugation.³⁶

Centrifugation should be carried out at ambient temperature (18-24°C). Temperature controlled centrifuges are not normally

required, but if used, temperatures of 2-8°C should not be used since this can affect results of some coagulation tests.¹²

Rapid deceleration at the end of centrifugation should be avoided when processing samples for coagulation testing to prevent remixing of cells and plasma around the plasma/buffy coat interface. Therefore, centrifuge braking should not be severe and use of swing out rotors is recommended to help prevent remixing.

If testing is performed on fresh plasma, most clot-based tests including PT, APTT and clotting factor assays can be performed on single centrifuged plasma. If samples are stored deep frozen prior to analysis, there are a number of tests where it is essential to remove as many platelets as practically possible prior to freezing, since freeze and thaw ruptures platelet membranes. Such rupture leads to release of platelet constituents, including platelet factor 4 (PF4) which can neutralize UFH in the sample and give rise to false low APTT or anti-Xa levels when monitoring heparinized patients.⁷ In addition, platelet membrane phospholipids are externalized with potential false normalization of some lupus anticoagulant (LA) tests.³⁷ Some tests for activated protein C resistance (APC-R) may also be affected by freezing once-centrifuged plasma.³⁸

To ensure that samples are sufficiently platelet-poor, they should be double centrifuged ("double-spun"), prior to freezing. To double spin plasma, the plasma from the first spin (not to include any of the cellular pellet if present) is aliquoted into a suitable noncontact activating secondary container (such as polypropylene) and re-centrifuged. The plasma from the second spin (not to include any of the cellular pellet if present) is aliquoted into an appropriate polypropylene tube and frozen.

Micropore filters that selectively remove some of the FVIII and VWF and other clotting factors from plasma should not be used to filter plasma in an effort to remove platelets.³⁹

Recommendation 5.1: Citrated blood samples for coagulation testing should be centrifuged at ambient temperature of 18-24°C using conditions sufficient to achieve platelet-poor plasma with a residual platelet count of $<10 \times 10^9/L$.

Remark: This can be normally be achieved by centrifugation for 15 minutes at 1500 g or 10 minutes at 1700 g and should be checked by performing platelet counts on platelet-poor plasma every 6-12 months.

Recommendation 5.2: Samples destined for lupus anticoagulant testing and samples containing UFH destined for APTT or anti-Xa testing should be subjected to a second centrifugation process after the first, prior to being frozen.

6 | PRIMARY TUBE STORAGE AND STABILITY

Clotting factors are potentially labile depending on the conditions of storage, so the interval between collection and processing as well as the storage conditions can affect the results of coagulation tests. There are many studies evaluating the stability of whole blood samples and centrifuged plasma for tests of haemostasis.⁴⁰⁻⁴⁸

It may be unsafe to generalize the findings from one study to other settings because the results may be different when different blood sample collection tubes are used, or when the tests are performed using different methods. The 5th edition of Clinical Laboratory Standards Institute (CLSI) guidance document H21 (under revision at the time of writing) recommends that most tests of haemostasis are performed within 4 hours of blood collection.⁸ There are, however, other publications addressing coagulation practice identifying tests/collection tubes/conditions associated with stability longer than 4 hours (see below). Where data are available, confirming additional stability for the specific tube and method in use, it is acceptable to extend the interval between collection and processing beyond 4 hours subject to any local regulatory requirements.

The stability of coagulation samples depends on a number of variables, including the blood collection system, whether the samples are processed prior to storage, storage temperature and the test parameters to be analysed.⁸ In general, whole blood samples should be kept capped and maintained at ambient temperature of 18-24°C prior to processing.

Some studies have demonstrated that citrated whole blood samples stored at ambient temperature are stable (defined by the article authors as <10% change) for up to 24 or 48 hours for many haemostasis tests including FII, FVII, FIX, FX, FXI, VWF:RCo, VWF:Ag, AT, PC, APC-R and D-Dimer (although not for APTT, factor V activity, factor VIII activity, and protein S activity).⁴⁷ Other studies have described significant changes in some test results over such time periods.⁴² Further evidence of D-Dimer stability for 24 hours in whole blood is available.⁴⁹ In some studies, citrated samples were stable for PT/International Normalized Ratio determinations at room temperature for 24 hours⁵⁰ or even longer.⁴⁰ Routine APTT determinations should be performed within 4 hours, although local studies may allow these limits to be exceeded.^{42,47} If testing is not to be performed within 4 hours for the APTT and 24 hours for PT, the plasma should be separated from the cellular fraction. Once separated, plasma can generally be maintained at room temperature or refrigerated for a few hours without adverse effects on results.⁵⁰ Factors II, VII, IX, X and FXI were found to be stable for at least 24 hours (defined as <10% change from baseline) in double centrifuged citrated plasma after transfer into polypropylene tubes with up to 12% reduction in FV and up to 14.8% reduction in FVIII after 4 hours.⁴⁸

Protein S activity is stable for 8 hours, whereas protein C, anti-thrombin (AT) and vitamin K dependent pro-coagulant factors are stable for at least 24 hours in plasma.^{41,42}

Samples containing UFH destined for APTT or anti-Xa testing must be centrifuged within 1 hour of collection and low molecular weight heparin containing samples within 24 hours of collection.^{50,51} PF4 released from platelets during storage of citrated samples will neutralize UFH, while LMWH is less susceptible to PF4 neutralization. Once centrifuged, APTT or anti-Xa testing on samples containing UFH should be completed within 4 hours of sample collection.⁵²

Sometimes, laboratories report the initially requested coagulation results and are then contacted later with a request to perform

an additional test on the original sample. It has been demonstrated that there was less than 5% change in PT, AT, fibrinogen or D-Dimer in plasma which remained over the cells in the primary tube after centrifugation for 24 hours at room temperature, although APTT may increase by close to 10%.⁵³

Recommendation 6.1: The allowed interval between sample collection and testing should be validated or verified for blood collection tube and specific coagulation laboratory tests and methods in use.

Recommendation 6.2: Citrated blood samples containing unfractionated heparin destined for APTT or anti-Xa determination should be centrifuged within 1 hour of collection, and the analysis shall be completed within 4 hours of collection.

7 | EFFECTS OF INTERFERING SUBSTANCES INCLUDING HAEMOLYSIS, ICTERUS AND LIPAEMIA

7.1 | Haemolysis

Visible haemolysis in the plasma of a centrifuged citrated blood sample may be caused by in vivo red cell disruption as a consequence of a particular disease process (eg, Auto-immune haemolytic anaemia, sickle cell crisis, transfusion reaction etc). It can also occur during extracorporeal circulation of the blood, or it may occur following in vitro lysis as an artefact of sample collection transport, storage or processing. If haemolysis is thought to have occurred in vivo, it is appropriate to proceed with analysis and reporting of coagulation tests. Results obtained on samples where haemolysis has occurred in vitro may be altered to the extent that they are unsafe to be used for patient management decisions (as described below), so that it is important to have data assessing the impact of haemolysis on coagulation test results.

It is usually considered that haemolysis is visible for plasma haemoglobin concentrations between 0.2 and 0.3 g/L or more,⁵⁴ although analytical systems used for haemoglobin measurements on such studies are sometimes designed for use with whole blood samples. Therefore, the recommended measuring range of the method may not include such low levels. Visual scrutiny of samples may be operator-dependent, and hence, it is useful that automated detection of haemolysis is available on multiple coagulation analysers.⁵⁵⁻⁵⁹ Automated assessment should be used where possible to standardize sample checks for presence of in vitro haemolysis. The optical characteristics and colour of plasma may occasionally be altered by therapeutic infusions which could impair interpretation of both visual and automated checks.

It is possible to perform a check for haemolysis in a blood sample in a few seconds using a point of care device which would allow the possibility of checking a tube while the collection device remains in the patient's vein so that a replacement tube could be collected before conclusion of the collection procedure without requirement for a new venipuncture.⁶⁰ This system could also be used when sample

analysis is being performed on an analyser which lacks automatic checks for haemolysis.

It has been common practice to assess the effects of haemolysis on coagulation tests by adding haemolysate to plasma, although it should be kept in mind that adding lysed red cells will not normally fully mimic all the effects of haemolysis when haemolysis occurs during blood sample collection and processing. There are a number of techniques in use to create spurious haemolysis including spiking plasma with haemolysate preparations; freezing and thawing whole anticoagulated blood; lysis of whole anticoagulated samples by means of deionized water with or without detergents; mechanical lysis of whole anticoagulated blood by sonication; stirring with a metallic bar; application of the blade of a tissue homogenizer and aspiration through a fine blood collection needle (<25 gauge).

There are a number of published studies describing the effects of haemolysis on results of clotting tests, in which cells have been deliberately lysed followed by addition of cell lysate to plasma. Addition of cell lysate to test plasma leads to dose-dependent trends to over-estimate PT and D-dimer with decrease in APTT and fibrinogen compared to baseline sample with no lysate added.⁶¹ This study reported a 5% increase in D-dimer results in the presence of 0.9% cell lysis, corresponding to 1.7 g/L plasma haemoglobin, with an increase of around 10% in D-dimer results at around 9 g/L plasma haemoglobin. For all the tests studied, there was considerable and unpredictable interindividual variation in the effects caused by addition of cell lysates to plasma prepared from the same normal subject. Significant differences in results of clotting tests as a consequence of targeted mechanical haemolysis have also been reported in other studies.^{55,62}

Interference in results of coagulation tests can derive from absorbance of the interfering substance at wavelengths that are commonly used for optical detection in coagulation analysers.⁶³ Beside this analytical interference, the interfering substance or interfering substance-related components can induce a direct interference with haemostasis. For instance, phospholipid membranes from haemolytic red cells may interfere with the coagulation reaction by providing a phospholipid-rich surface that accelerates coagulation reactions.⁶⁴ Platelet activation by haemolysed red cells can impact coagulation test results indicating that effects will not be restricted to those induced by change in the colour of the plasma samples under analysis. The disruption of cells may not be restricted to erythrocytes and release of cellular content will not be restricted to haemoglobin. It has been shown that red blood cell lysis activates platelets through ADP release and that the haemoglobin released enhances platelet activation by lowering nitric oxide bioavailability.²⁵ These effects are likely to contribute to changes in coagulation parameters in haemolysed blood samples during the interval between occurrence of haemolysis and analysis of the sample in the laboratory.

Haemolysis in citrated blood samples which occurs *in vitro* is normally preventable if appropriate collection transport and processing techniques are used, but nevertheless is commonly encountered. It is unclear to what extent the impact of haemolysate added to plasma reproduces all the effects on coagulation tests that may

TABLE 1 Causes of *in vitro* haemolysis (adapted from reference 63)

Sample collection	Poor venous access
	Poor technique
	Traumatic venipuncture
	Small needle size
	Forced passage through needle
Sample handling and transportation	Excessive shaking or trauma to tube
	Excessive delay before centrifugation
	Exposure to very high or very low temperatures
Sample preparation	Inappropriate centrifugation conditions, for example temperature too high or low.
	Transfer of plasma contaminated by red blood cells followed by freeze/thaw

occur in samples where haemolysis has occurred during sample collection and processing, particularly because there are a number of causes of unintended haemolysis in citrated blood samples which have been reviewed elsewhere⁶³ and are summarized in Table 1.

Since there are many possible causes of *in vitro* haemolysis, it is likely that the effects on clotting tests are also multifactorial with coagulation activation, clotting factor consumption and changes in the optical and mechanical characteristics of the sample occurring in isolation or in combination in different samples. Some of these effects may be absent in studies of the effect of haemolysis on coagulation tests performed using samples which have been constructed by addition of lysate to plasma. Such studies should be interpreted with caution.

Some centres reject samples with *in vitro* haemolysis and request replacements. Some studies have therefore been able to compare results in rejected haemolysed samples with results obtained on subsequently received replacement nonhaemolysed samples from the same subjects collected within 3 hours⁶² or within 4 hours of each other.⁵⁹

The level of free plasma haemoglobin in samples with visible *in vitro* haemolysis was 0.2-7.0 g/L in one of these and 0.5-9 g/L in another (total sample number in the 2 studies was n = 119). The two studies used different blood collection tubes from widely used international suppliers, suggesting that the level of *in vitro* haemolysis is not a consequence of poor collection tube quality and not specific to a particular tube type.

Both studies show that, on average, the impact of haemolysis on PT and APTT was low but there were individual samples in which changes in APTT induced by *in vitro* haemolysis would have important clinical implications. This did not occur for PT, that is the bias in haemolysed samples were not clinically relevant. Falsely normal APTTs occurred in approximately 5% of haemolysed samples where the true APTT (on a matched nonhaemolysed sample) was more than 6 seconds prolonged but was normal on the haemolysed sample.⁶⁵ In other samples, D-Dimer results that should have been below the cut-off for venous thromboembolism exclusion were falsely elevated

above the cut-off in haemolysed samples.⁶⁵ Detection and rejection of haemolysed samples for APTT and D-Dimer is therefore important for patient management purposes. Clinically relevant differences in APTT and D-Dimer, but not PT or Antithrombin, were reported in patient samples with spontaneous haemolysis when compared to matched nonhaemolysed samples from the same subjects using a different photo-optical analytical platform to the studies described above.⁵⁸ If a request for D-Dimer is made to confirm a level is below the cut-off for VTE, then it may be possible to safely test a haemolysed sample, despite some possible artificial elevation, provided the result remains below the cut-off, since D-Dimer is not falsely reduced in these studies (being either unchanged or falsely elevated). Here, a nonclinically significant (artificial) increase caused by haemolysis may have no impact on management whereas rejection could introduce delay.

The APTT of a sample with in vitro haemolysis is often falsely shortened, but can also be falsely prolonged, and the change in APTT as a consequence of in vitro haemolysis is not correlated with the level of plasma haemoglobin irrespective of whether the APTT was determined photo-optically⁶³ or mechanically.⁵⁹ Clinically relevant changes in APTT can occur when even very low levels of in vitro haemolysis are present.^{59,64,65} It is therefore unsafe to attempt interpretation of APTT on samples with in vitro haemolysis.

In the studies comparing results obtained on matched pairs of haemolysed and nonhaemolysed samples from the same patient, the authors attributed the differences observed to the effects of haemolysis, but the possibility of a real change in the patients clotting results within the period between collection of the two samples (up to 4 hours in the studies discussed) could not be fully excluded in all cases. Despite this limitation, conclusions on the impact of in vitro haemolysis on tests of coagulation may be safer when based on data comparing haemolysed and nonhaemolysed samples from the same patient where possible.

Haemolysis frequently occurs in vivo in thrombotic microangiopathies, but if haemolysis is thought to have occurred in vitro as an artefact, in citrated samples for ADAMTS 13 activity determination, then ICSH have recommended that replacement samples should be obtained, and that if this is not possible, a comment regarding potential assay interference should be added to the laboratory report.⁶⁶

Recommendation 7.1.1: Plasma for coagulation testing which has been prepared from citrated blood samples should be checked for the presence of in vitro haemolysis, preferably using an automated system for consistency.

Recommendation 7.1.2: APTT should not be performed on samples with haemolysis that has occurred in vitro during sample collection, transport and processing.

Remark: The impact of moderate in vitro haemolysis (ie, <10 g/L) on PT/INR is usually clinically irrelevant.

Recommendation 7.1.3: The possibility that haemolysis has occurred in vivo should be considered. Samples from patients with in vivo haemolysis for determination of coagulation tests can be accepted and tested.

7.2 | Icterus

Icterus is the yellow colour caused by bilirubinaemia in plasma. Bilirubin has a high absorbance between 400 and 520 nm with a peak at 456 nm which has the potential to interfere when photo-optical coagulation analysers use similar wavelengths to monitor reactions.

There are a number of studies assessing interference of bilirubin on coagulation testing using photo-optical analysers, sometimes by comparing results before and after addition of bilirubin to plasma, and sometimes by comparison of results in patient samples with hyperbilirubinaemia obtained with optical and mechanical coagulometers including samples with markedly elevated concentrations (for a comprehensive review see 63).

The impact of icterus on PT, APTT or fibrinogen with a mechanical system was clinically irrelevant in one study,⁵⁹ while another study reported some shortening of PT and APTT (all <10%), increase in Fibrinogen (up to 20%) and elevation of D-dimer (<10%) in the presence of increasing bilirubin concentrations, depending upon reagent/instrument combination.⁶⁷ Bilirubin levels >100 µmol/L can reduce ADAMTS 13 activity in some fluorescence assays by quenching fluorescence which can sometimes be avoided by sample dilution.⁶⁶ Sample dilution may be an option in some other tests such as chromogenic assays, if bilirubin interference is suspected. Where a manufacturer identifies a level of bilirubin that causes interference in a particular test that should be taken into account when interpreting and reporting test results.

Recommendation 7.2.1: Citrated samples with hyperbilirubinaemia for 387726645PT APTT and fibrinogen coagulation tests should be accepted for analysis.

7.3 | Lipaemia

Lipaemia in citrated plasma samples can lead to problems when samples are analysed for coagulation tests due to increased light scatter or absorbance using photo-optical systems.⁶³ Lipid can be added to plasma samples to simulate the lipaemia that occurs in patient samples to assess the level of interference but it is not known to what extent these samples truly mimic patient samples where different types and concentrations of lipid may be present. The interference of lipid can be reduced by using increased centrifugal force.^{63,68}

One study compared results on samples before and after use of a double high speed centrifugation (often referred to as ultracentrifugation) of 15 000 g for 15 minutes in 101 samples with excess lipaemia (detected by an analyser). The differences in mean results of PT, APTT AT and D-Dimer were <4%. Extra centrifugation increased the result of Clauss fibrinogen measurements by approximately 10% for fibrinogen in the range 2-4 g/L, suggesting that there was some underestimation in lipaemic samples.⁶⁹ In another study, centrifugation at 10 000 g for 10 minutes was sufficient to markedly reduce the lipid to the extent that all 22 samples failing to give PT or APTT results in the presence of lipaemia (mix of constructed and real patient samples) gave clotting times after centrifugation.⁷⁰

The impact of lipaemia on PT, APTT or fibrinogen with a mechanical system was clinically irrelevant in one study⁵⁹ and dependent upon the reagent/instrument combination in another study.⁶⁷

Samples with detectable lipaemia should not be rejected since replacement samples may have similar characteristics and because the impact of excess lipid found in patient samples on coagulation test results is small and not normally clinically relevant.

Recommendation 7.3.1: If the lipid level is so high that the analyser cannot accurately detect clotting, the sample can be analysed after removal of the lipid from the plasma by centrifugation at 10 000 g for 10 minutes at ambient temperature.

Remark: If high speed centrifugation is used, a reference range should be established for tests on samples processed using that procedure.

7.4 | Other interfering substances

It has been reported that C-reactive protein (CRP) can cause prolongation of APTT depending on the reagent used. There was a positive correlation between increasing CRP concentration and APTT determined with STA Cephascreen APTT reagent but not with a kaolin-based APTT reagent in 59 emergency department patients. Addition of CRP to normal plasma caused average increases in APTT with four reagents as follows: STA Cephascreen 5.6 seconds; STA CK Prest 1.6 seconds; Actin FS 3.5 seconds; Synthasil 2.1 seconds.⁷¹ There is also some evidence that CRP can interfere in lupus anticoagulant testing.⁷² A number of other substances may interfere in coagulation tests. Interfering substances should be considered when unexpected results occur and repeat sample collection may be useful in such case.

8 | SECONDARY ALIQUOTS—TRANSPORT, STORAGE AND PROCESSING

When samples are submitted to the laboratory as secondary tubes, samples collected into the wrong anticoagulant or that have been combined from blood collected into different types of tube are not readily recognized as inappropriate samples. To avoid reporting results of tests performed on an inappropriate matrix, some laboratories characterize (by measuring different electrolytes) those samples that yield abnormal results and have been submitted to the laboratory as secondary tubes before reporting.⁷³ See ICSH recommendations for collection and processing of blood samples for coagulation testing for a review of coagulation test results on samples collected as serum or into anticoagulants other than citrate.⁷

If tests for lupus anticoagulant or UFH (APTT or anti-Xa for example) are to be performed on plasmas which have been frozen prior to analysis, the plasma should be separated into a secondary tube after the first centrifugation and subjected to a second centrifugation process, as described in Section 5, prior to freezing.

Secondary tubes constructed from polypropylene are suitable for storage of frozen citrated plasma for coagulation test.⁷⁴ Differences

between results obtained on samples stored in screw cap tubes with rubber seals and snap top versions were clinically irrelevant and either can be used.⁷⁴ Secondary aliquots for deep freezing should be capped throughout preparation, freezing and thawing.⁷⁴

When platelet-poor plasma (PPP) is frozen for future testing, the freezing and storage conditions affect the stability of frozen material. Rapid freezing reduces the opportunity for changes to the clotting factors and is more easily achieved when the secondary tube is transferred directly into storage units already present in the freezer for sufficient time to be acclimatized to the freezer temperature, rather than transferring samples in carriers or storage units that have been at ambient temperature.

Storage of platelet-poor plasma (PPP) at -20°C may be inadequate to prevent change in results of some coagulation tests. Freezers with autodefrost should not be used to store PPP prior to determination of coagulation tests since there can be significant rise in temperature during defrost.

The sample stability is improved at lower temperatures and while -24°C is adequate for 3 months storage, stability for many coagulation analytes is improved by use of -70°C for samples that are being stored for more than 3 months.⁷⁵ Minor differences between results on samples tested within 2-4 hours of collection and again after 1 week storage at -70°C were not clinically significant for PT APTT, fibrinogen, D-Dimer, AT, FV, FVII, FVIII and FIX.⁷⁶

If frozen plasma samples are transferred to another site for analysis, this should be done using conditions which maintain the sample temperature below -20°C. This is most effectively done by packaging the samples in close contact with sufficient volume of solid carbon dioxide. Use of temperature loggers alongside such samples provides useful evidence that adequately low temperatures have been maintained. Secondary samples transported like this must utilize sample containers with adequate sealing of the cap since CO₂ gas entering the samples can lead to pH change in lyophilized plasmas, to the extent that some clotting test results could be affected.⁷⁷ Other studies have reported pH changes in plasma samples exposed to dry ice.^{78,79}

Prolongations of PT and APTT in frozen citrated plasma which occurred following 24 hours exposure to dry ice were abolished if such plasmas were stored at -80°C for 24 hours between dry ice exposure and analysis.⁷⁸

Frozen citrated plasma samples should not be allowed to thaw slowly at room temperature or precipitation of fibrinogen, FVIII, VWF and AT may occur. Samples should be thawed in a temperature controlled water bath at 37°C with the surface of the frozen plasma held at or below the surface of the water, after which the sample must be mixed by 3-5 over end inversions prior to testing. If the volume of plasma is less than 1 mL, then 5 minutes thawing time is sufficient.⁷⁴ Thawing appropriately (at 37°C for 5 minutes) without post-thaw mixing is not sufficient since a lack of mixing was associated with a 60% loss of VWF activity and antigen from normal plasmas which was avoided by any of 3 different mixing methods, that is gentle mixing by 6-8 inversions; 5 minutes on a roller mixer or vortexing for 10 seconds.⁸⁰

It has historically been common practice to avoid refreezing of secondary plasma aliquots which have been previously frozen and

then thawed despite a lack of published studies supporting that practice. Actually, there are several publications addressing this issue, which indicate that some coagulation parameters are unaffected by multiple freeze and thaw cycles, at least under the conditions studied. Using 10% change as significant, FII, FV and FIX were not significantly affected by two freeze-thaw cycles whereas FXII was significantly affected by a second freeze thaw in one study.⁸¹ A comprehensive study assessing seven cycles of freeze thaw (with -70°C storage) confirmed less than 5% change after two cycles for all parameters studied and <6% change after 3 cycles.⁷⁴ The authors studied PT, APTT, Fibrinogen, D-Dimer, clot-based FII/FV/FVII/FVIII/FIX/FXI/FXII, chromogenic FXIII, AT activity, PC activity, PS activity plasminogen, DRVVT and anti-Xa (LMWH). Five freeze-thaw cycles did not cause significant changes in anticardiolipin or anti β 2 glycoprotein 1 antibody testing.⁸² Taken together, these studies indicate that many coagulation tests can be safely performed after 2 cycles of freeze thaw when thawing is done in a water bath for 4-5 minutes and freezing/refreezing is done at -70°C, provided the interval between thaw and refreeze takes account of sample stability (see Section 6 above).

Recommendation 8.1 Citrated plasma should be frozen in polypropylene vials which are capped at all stages of processing, freezing and thawing.

Recommendation 8.2: Citrated platelet-poor plasma samples can be stored at -24°C for up to 3 months or -70°C for at least six months prior to performing coagulation tests.

Recommendation 8.3: If frozen citrated plasma samples are transported between sites, they should be maintained at -20°C or colder during transfer.

Recommendation 8.4: Frozen citrated plasma samples of up to 1 mL in volume must be thawed rapidly in a water bath for 4-5 minutes at 37°C and thoroughly mixed to avoid formation of cryoprecipitate.

Recommendation 8.5: Use of 2 or more cycles of freeze thaw prior to performing coagulation tests may be acceptable, but should be verified for the local samples and methods being used.

9 | PREANALYTICAL VARIABLES FOR PLATELET FUNCTION TESTING

There are several guidelines, consensus documents and review articles that provide useful information about platelet function testing (PFT).⁸³⁻⁸⁵ Preanalytical variables that impact on PFT in the assessment of aspirin and clopidogrel responsiveness have been reviewed⁸⁶ and will not be specifically discussed here.

9.1 | Blood collection for platelet function testing

The collection and handling of blood samples for PFT requires special attention in order to reduce platelet activation during handling. The specific details may vary depending on the reason for the study

(investigation of inherited and acquired platelet function defects, assessment of responsiveness to platelet antagonist treatment) and the type of platelet function test to be performed (eg, PFA-200, light transmission aggregometry [LTA], impedance aggregation [IA], flow cytometric techniques, thromboelastography).

Where possible, a rest period of few minutes before blood collection may reduce stress and exercise induced changes. A record of all drugs (including nonprescription or commonly referred as over-the-counter medicines) and herbal supplements that have been taken in the last 7 days should be collected, due to the large number of substances that may affect PFT results (see 86 for a useful list).

For most purposes, it is generally recommended that blood for platelet function is collected by clean venipuncture, using a needle size of at least 21 gauge, rapidly anticoagulated (usually in citrate) using siliconized glass or polypropylene blood collection tubes, maintained at room temperature (18-24°C; never chilled) and tests completed within 4 hours. National and International consensus and guidelines have recommended either 109 or 129 mmol/L buffered sodium citrate for LTA,⁸⁵⁻⁸⁷ but there is some evidence that 109 mmol/L is preferable and should be tested within 2 hours.^{88,89} It is also suggested that the first 3-4 mL of blood is discarded or used for other purposes.⁷ There is some evidence that citrate tube volume may also have an effect, with a decrease in impedance aggregation with time, using collection tubes less than 3.6 mL.⁹⁰

For some IA methods, hirudin and BAPA anticoagulated blood have been recommended, as these may improve the storage stability of the sample in relation to certain agonists.⁹¹⁻⁹⁴ However, platelet clumping may be observed in hirudin anticoagulated blood with a consequent decrease in IA at 2-3 hours.⁹⁵ The choice of anticoagulant, specimen collection, transport and time between sampling and analysis all have critical effects on the closure time using the PFA-series analysers.⁹⁶⁻⁹⁹

For flow cytometric investigations of platelet membrane glycoprotein receptors or generation of activation markers in response to agonists, citrated whole blood is normally used. A recent study showed no differences between vacuum tube and syringe aspirated samples or between the first and second tubes drawn.¹⁰⁰ Fluorophore-labelled antibodies should be added to the sample as soon as possible after blood collection; after appropriate incubation and dilution, a mild fixative should be added to prevent further activation or deterioration (for guidance and protocols see ¹⁰¹⁻¹⁰³).

The impact of blood collection method has recently been studied and found to have no significant effect on PFA and IA (using the Multiplate analyser) methods. Blood collection can be by venipuncture and from venous, arterial and central venous lines.¹⁰⁴ However, this may depend on the type of substances that have previously been infused through the line (ie, blood products or heparin) and there may be differences between research studies with carefully collected samples and the real world, with blood collection outside controlled conditions. Some of the authors also found no difference in PFA results for syringe collection versus evacuated tube.

9.2 | Transport of samples for platelet function testing

Transporting blood samples to the laboratory by pneumatic tube systems (PTS) may considerably reduce the turnaround time for testing, but each type of system should be assessed for impact on PFT results. In Section 3 above, we have recommended against use of PTS to transport blood samples prior to tests of platelet function although this may be possible in specific circumstances if robust data are available to support the use of particular defined conditions since one recent study found no difference in PFT results (assessed by LTA and flow cytometric methods) when a PTS was compared to the regular transport system.¹⁰⁵ However, others have found significant changes in LTA and pre-activation of platelets by flow cytometry¹⁰⁶ decreased platelet function by PFA and IA methods.^{107,108} See also Section 3 above.

9.3 | Sample preparation for platelet function testing

Blood samples should be allowed to “rest” at room temperature for 15 minutes before centrifugation for platelet rich plasma (PRP). Some centres utilize a similar “rest” period for PRP before testing, in order to allow restabilization after the mechanical trauma of centrifugation. PRP should be prepared by centrifugation at 200-250 g for 10 minutes at room temperature, as loss of function occurs in LTA with increasing centrifugal force, probably due to the loss of larger platelets as indicated by a falling platelet count and decreasing mean platelet volume.^{109,110}

Haemolysed samples should not be used (unless there is good reason to expect in vivo haemolysis), since this implies that the sample is activated and therefore platelet function will be altered in vitro.²⁵ Samples containing small clots or clumped platelets should not be tested. Grossly lipaemic samples may cause problems in some optical methods (eg, LTA, where it may be difficult to set a suitable 100% light transmittance).

There is controversy about the need to adjust the PRP platelet count for LTA. The general consensus is that adjustment is not necessary when the PRP count is within the normal range as this can impair responsiveness to agonists.^{85,111} However, most workers feel that the platelet count should be reduced to normal when above $600 \times 10^9/L$ and this is usually performed using autologous platelet-poor plasma (PPP; prepared by centrifuging at 2000 g for 15 minutes). PRP and PPP should be handled using polypropylene pipettes and tubes and should be stored at room temperature in capped tubes to reduce pH changes.

Recommendation 9.1: Blood for platelet function testing should be collected using a 19-21 gauge needle, mixed with 109 mmol/L trisodium citrate (excluding the first 3-4 mL collected) and maintained at 18-24°C during processing and testing.

Recommendation 9.2: Platelet rich plasma should be prepared by centrifugation for 10 minutes at 200-250 g at room temperature of 18-24°C. The platelet count of PRP should only be adjusted down

to $250-300 \times 10^9/L$ when the platelet count exceeds $600 \times 10^9/L$. Any such adjustment should use autologous platelet-poor plasma.

Recommendation 9.3: Processing of blood samples should begin 15-30 minutes after collection and analysis must be completed preferably within 2 hours and always within a maximum of 4 hours, taking account of any relevant manufacturers' instructions for use.

10 | CONCLUSIONS

Some blood samples may not be easy to replace, and not all rejected samples are replaced within the relevant time frame. This means that any citrated blood sample acceptance/rejection protocols taking account of the causes of bias in test results described above and in ICSH recommendations on sample collection⁷ should balance the risks of releasing unsafe or misleading results against the risks associated with rejection followed by delay or absence of coagulation test results. If a preanalytical change causes bias, which is unequivocally not clinically relevant, then the balance of risk favours release of those results. This type of approach has been considered in other areas of laboratory medicine.¹¹²

Recommendation 10.1: All coagulation laboratories should establish written policies on sample transport, storage and processing of both primary tubes and secondary aliquots which includes a description of which citrated samples are acceptable for testing and those that must be rejected; this should be jointly developed by laboratory management and the clinicians who make use of the coagulation laboratory service.

ACKNOWLEDGEMENTS

The drafting group are grateful to the Western Norway Regional Health Authority for supporting Ann Helen Kristoffersen with post-doctoral fellowships to facilitate participation in drafting this guidance.

CONFLICT OF INTEREST

The authors have no competing interests.

DATA AVAILABILITY STATEMENT

Data sharing is not applicable to this article as no new data were created or analyzed in this study.

ORCID

Steve Kitchen  <https://orcid.org/0000-0003-4790-1159>
Dorothy M. Adcock  <https://orcid.org/0000-0002-0182-6272>
Giuseppe Lippi  <https://orcid.org/0000-0001-9523-9054>
Ian Mackie  <https://orcid.org/0000-0002-7393-2365>
Richard A. Marlar  <https://orcid.org/0000-0002-8223-2811>
Sukesh Nair  <https://orcid.org/0000-0002-3760-0929>

REFERENCES


1. Lippi G, Franchini M, Montagnana M, Salvagno GL, Poli G, Gudi GC. Quality and reliability of routine coagulation testing: can we trust that sample? *Blood Coagul Fibrinolysis*. 2006;17:513-519.

2. Favaloro EJ, Lippi G, Adcock DM. Pre analytical and post analytical variables: the leading causes of diagnostic errors in haemostasis? *Semin Thromb Hemost*. 2008;34:612-634.
3. Kalra J. Medical errors: impact on clinical laboratories and other clinical areas. *Clin Biochem*. 2004;37:1052-1062.
4. Bonini P, Plebani M, Ceriotti F, Rubboli F. Errors in laboratory medicine. *Clin Chem*. 2002;48:691-698.
5. Green SF. The cost of poor blood specimen quality and errors in preanalytical processes. *Clin Biochem*. 2013;46:1175-1179.
6. Salvagno GL, Lippi G, Bassi A, Poli G, Guidi GC. Prevalence and type of pre analytical problems for in patient samples in coagulation laboratories. *J Eval Clin Pract*. 2008;14:351-353.
7. Kitchen S, Adcock DM, Dauer R, et al. International Council for Standardisation in Haematology (ICSH) recommendations for collection of blood samples for coagulation testing. *Int J Lab Hematol*. 2021;43(4):571-580. 10.1111/ijlh.13584
8. Clinical Laboratory Standards Institute. *Collection, Transport and Processing of Blood Specimens for Testing Plasma-Based Coagulation Assays And Molecular Haemostasis Assays: Approval Guidelines*, 5th edition. CLSI H21-A510; 2008.
9. Kim YA, Lewandrowski KB, Lucien FA, Van Cott EM. The effects of transport temperature and time on routine and specialized coagulation assays. *Blood Coagul Fibrinolysis*. 2018;29:184-188.
10. Favaloro E, Nair SC, Forsyth CJ. Collection and transport of samples for laboratory testing in von Willebrand's disease (VWD): time for a reappraisal? *Thromb Haem*. 2001;86:1589-1590.
11. Kitchen S, Malia RG, Greaves M, Preston FE. A method for the determination of activated factor VII using bovine and rabbit thromboplastins: demonstration of increased levels in disseminated intravascular coagulation. *Thromb Res*. 1988;50:191-200.
12. Favaloro EJ, Soltani S, McDonald J. Potential laboratory misdiagnosis of hemophilia and von Willebrand disorder owing to cold activation of blood samples for testing. *Am J Clin Pathol*. 2004;122:686-692.
13. Böhm M, Täschner S, Kretzschmar E, Gerlach R, Favaloro EJ, Scharer I. Cold storage of citrated whole blood induces drastic time-dependent losses of factor VIII and von Willebrand factor: potential for misdiagnosis of haemophilia and von Willebrand disease. *Blood Coagul Fibrinolysis*. 2006;17:39-45.
14. Refaii MA, van Cott EM, Lukoszyk M, Hughes J, Eby CS. Loss of factor VIII and von Willebrand activities during cold storage of whole blood is reversed by rewarming. *Lab Hematol*. 2006;12:99-102.
15. Amukele TK, Sokoll LJ, Pepper D, Howard DP, Street J. Can unmanned aerial systems (Drones) be used for routine transport of chemistry, haematology and coagulation laboratory specimens? *PLoS One*. 2015;10(7):e0134020. 10.1371/journal.pone.0134020
16. Kratz A, Salem RO, Van Cott EM. Effects of a pneumatic tube system on routine and novel haematology and coagulation parameters in healthy volunteers. *Arch Pathol Lab Med*. 2007;131:293-296.
17. Dyszkiewicz-Korpanty A, Quinton R, Jassine J, Sarode R. The effect of pneumatic tube transport on PFA-100TM closure time and whole blood aggregation. *J Thromb Haemost*. 2004;2:354-356.
18. Heireman L, Van Geel P, Musger L, Heylen E, Uyttenbroeck W, Mahieu B. Causes, consequences and management of sample haemolysis in the clinical laboratory. *Clin Biochem*. 2017;18:1317-1322.
19. Gils C, Broell F, Vinholt PJ, Nielsn C, Nyobo M. Use of clinical data and acceleration profiles to validate pneumatic transportation systems. *Clin Chem Lab med*. 2020;58:560-568.
20. Guss DA, Chan TC, Killeen JP. The impact of a pneumatic tube and computerised physician order management on laboratory turnaround time. *Ann Emerg Med*. 2008;51:181-185.
21. Nybo M, Lund ME, Tittlestad K, Maegaard CU. Blood sample transportation by pneumatic transportation systems: a systematic literature review. *Clin Chem*. 2018;64:782-790.
22. Pragay DA, Edwards L, Toppin M, Palmer RR, Chilcote ME. Evaluation of an improved pneumatic tube system suitable for transportation of blood specimens. *Clin Chem*. 1974;20:57-60.
23. Cui M, Jing R, Wang H. Changes of serum lactate dehydrogenase and potassium levels produced by a pneumatic tube system. *Lab Med*. 2009;40:728-731.
24. Gomez-Rioja R, Fernandez-Calle P, Alcaide MJ, et al. Inter-individual variability of hemolysis in plasma samples during pneumatic tube system transport. *Clin Chem Lab Med*. 2013;51:e231-e233.
25. Helms CC, Marvel M, Zhao W, et al. Mechanisms of hemolysis-associated platelet activation. *J Thromb Haemost*. 2013;11:2148-2154.
26. Streichert T, Otto B, Schnabel C, et al. Determination of hemolysis thresholds by the use of data loggers in pneumatic tube systems. *Clin Chem*. 2011;57:1390-1397.
27. ISO 15189. Medical laboratories – Requirements for quality and competence; 2012.
28. Farnsworth CW, Webber DM, Krekler JA, Budelier MM, Bartlett NL, Gronowski AM. Parameters for validating a hospital pneumatic tube system. *Clin Chem*. 2019;65:694-702.
29. Cooper J, der Van Hoop J, Broell F, Kitchen S. Identifying and reducing high levels of haemolysis from samples transport in pneumatic tube systems. *J Thromb Haemost*. 2020. In press.
30. Delbruck C, Miesbach W. The course of von Willebrand factor and factor VIII activity in patients with von Willebrand disease during pregnancy. *Acta Haematol*. 2019;142:71-78.
31. Austin AW, Wirtz PH, Patterson SM, Stutz M, von Kanel R. Stress-induced alterations in coagulation: assessment of a new hemocoagulation correction technique. *Psychosom Med*. 2012;74:288-295.
32. Bowyer A, Kitchen S, Maclean R. Effects of emicizumab on APTT, one-stage and chromogenic assays of factor VIII in artificially spiked plasma and in samples from haemophilia A patients with inhibitors. *Haemophilia*. 2020;26:536-542.
33. Lippi G, Salvagno GL, Monyagna M, Manzaato F, Guidi GC. Influence of the centrifuge time of primary plasma tubes on routine coagulation testing. *Blood Coagul Fibrinol*. 2007;18:525-528.
34. Carroll WE, Wollitzer AO, Harris L, Ling MC, Whitaker WL, Jackson RD. The significance of platelet count in coagulation studies. *J Med*. 2001;32:83-86.
35. Nelson S, Pratt A, Marlar RA. Rapid preparation of plasma for 'stat' coagulation testing. *Arch Pathol Lab Med*. 1994;118:175-176.
36. Kao CH, Shy LC, Yen WF. Evaluation of a high-speed centrifugation with rapid preparation of plasma for coagulation testing to improve turnaround time. *J Biomed Lab Sci*. 2010;22:23-27.
37. Pengo V, Tripodi A, Reber G, et al. Subcommittee on lupus anti-coagulant/antiphospholipid antibody of the scientific and standardisation committee of the International Society on Thrombosis and Haemostasis. Update of the guidelines for lupus anticoagulant detection. *J Thromb Haemost*. 2009;7:1737-1740.
38. Tripodi A, Valsecchi C, Chantarangkul V, Battaglioli MPM. Standardisation of activated protein C resistance testing: effects of residual platelets in frozen plasmas assessed by commercial and home-made methods. *Br J Haem*. 2003;120:825-828.
39. Favaloro EJ. Pre analytical variables in coagulation testing. *Blood Coagul Fibrinolysis*. 2007;18:86-89.
40. Baglin T, Luddington R. Reliability of delayed INR determination: implications for decentralised anticoagulant care with off-site blood sampling. *Br J Haem*. 1997;96:431-434.
41. Heil W, Grunewald R, Amnd M, Heins M. Influence of time and temperature on coagulation analysis in stored plasma. *Clin Chem Lab Med*. 1998;36:459-452.
42. Awad MA, Selim TE, Al-Sabbagh FA. Influence of storage time and temperature on International Normalized Ratio (INR) levels

- and plasma activities of vitamin K dependent clotting factors. *Hematology*. 2004;9:333-337.
43. O'Neill EM, Rowley J, Hanson-Wicher H, McCarter S, Ragno G, Valeri CR. Effect of 24-hour whole-blood storage on plasma clotting factors. *Transfusion*. 1999;39:488-491.
 44. Rao LV, Okorodudu AO, Petersen JR, Elghetany MT. Stability of prothrombin time and activated partial thromboplastin time tests under different storage conditions. *Clin Chim Acta*. 2000;300:13-21.
 45. van Geest-Daaderop JH, Mulder AB, Boonman-deWinter LJ, Hoekstra MM, van den Besselaar AM. Preanalytical variables and off-site blood collection: influences on the results of the prothrombin time/international normalized ratio test and implications for monitoring oral anticoagulant therapy. *Clin Chem*. 2005;51:561-568.
 46. Toulon P, Metge S, Hangard M, Zwahlen S, Paiulenne S, Besson V. Impact of different storage times at room temperature of unspun citrated blood samples on routine coagulation test results. Results of a bicentre study and review of the literature. *Int J Lab Hematol*. 2017;39:458-468.
 47. Zurcher M, Sulzer I, Barizzi G, Lammle B, Alberio L. Stability of coagulation assays performed in plasma from citrated whole blood transported at ambient temperature. *Thromb Haemost*. 2008;99:416-426.
 48. Linskens EA, Devreese KMJ. Pre-analytical stability of coagulation parameters in plasma stored at room temperature. *Int J Lab Hematol*. 2018;40:292-303.
 49. Caliezi C, Reber G, Lämmle B, de Moerloose P, Wuillemin WA. Agreement of D-dimer results measured by a rapid ELISA (VIDAS) before and after storage during 24 hr or transportation of the original whole blood samples. *Thromb Haemost*. 2000;83:177-178.
 50. Adcock DA, Kressin DC, Marlar RA. The effect of time and temperature variables on routine coagulation tests. *Blood Coagul Fibrinolysis*. 1998;9:463-470.
 51. Birri N, Baumgartner D, Conte T, et al. Stability of low molecular weight heparin anti-factor Xa activity in citrated whole blood and plasma. *Br J Haematol*. 2011;155:620-641.
 52. Adcock Funk DM, Lippi G, Favaloro EJ. Quality Standards for sample processing transportation and storage in haemostasis testing. *Semin Thromb Hemost*. 2012;38:576-585.
 53. Rimac V, Coen HD. Is it acceptable to use coagulation plasma samples stored at room temperature and 40C for 24 hours for additional PT APTT fibrinogen, antithrombin and D-dimer testing? *Int J Lab Hematol*. 2017;39:475-481.
 54. Lippi G, Guidi GC, Mattiuzzi C, Plebani M. Pre-analytical variability: the dark side of the moon in laboratory testing. *Clin Chem Lab Med*. 2006;44:358-365.
 55. Tantanate C, Teyateeti M, Tiendதாகul P. Influence of plasma interferences on screening coagulation and performance evaluation of the automated coagulation analyser Sysmex CS2100. *Siriraj Med J*. 2010;63:151-156.
 56. Lippi G, Ippolito L, Favaloro EJ. Technical evaluation of the novel pre-analytical module on instrumentation laboratory ACL TOP: advancing automation in hemostasis testing. *J Lab Autom*. 2013;18:382-390.
 57. Florin L, Oyaert M, Van Maerken T, Devreese KMJ. Performance evaluation of the pre-analytical check module of the Stago STA R Max2 mechanical endpoint detection analyzer for the assessing of the impact of hemolysis lipemia and icterus on APTT and PT. *Int J Lab Hematol*. 2018;40:e109-e112.
 58. Montanuli B, Guitto C, Cosseddu D. Influence of hemolysis, icterus and lipemia on coagulation tests as performed on Cobas t511 new analyser. *Blood Coagul Fibrinol*. 2020;31:48-54.
 59. Woolley A, Golmard JL, Kitchen S. Effects of haemolysis, icterus and lipaemia on coagulation tests as performed on Stago STA-Compact-Max analyser. *Int J Lab Hematol*. 2016;38:375-388.
 60. Duhalde H, Skogo J, Karlsson M. Point-of-Care hemolysis detection in blood gas specimens directly at the emergency department. *Scan J Clin Lab Investig*. 2019;79:283-287.
 61. Lippi G, Montagnana M, Salvagno GL, Guidi GC. Interference of blood cell lysis on routine coagulation testing. *Arch Pathol Lab Med*. 2006;130:181-184.
 62. Lippi G, Avanzini P, Zobbi V, Ippolito L. Influence of mechanical hemolysis of blood on two D-dimer immunoassays. *Blood Coagul Fibrinolysis*. 2012;23:461-463.
 63. Lippi G, Plebani M, Favaloro EJ. Interference in coagulation testing: focus on spurious hemolysis, icterus, and lipemia. *Semin Thromb Hemost*. 2013;39:258-266.
 64. Laga AC, Cheves TA, Sweeney JD. The effect of specimen hemolysis on coagulation test results. *Am J Clin Pathol*. 2006;126:748-754.
 65. Woolley A, Horner K, Kitchen S. Haemolysis during sample collection causes variable but sometimes clinically relevant important errors. *J Thromb Haemost*. 2011;9(Suppl. 2):125.
 66. Mackie I, Mancini I, Muia J, et al. International Council for Standardization in Hematology (ICSH) recommendations for laboratory measurement of ADAMTS13. *Int J Lab Hematol*. 2020;42:685-696.
 67. Nagant C, Rozen L, Demulder A. HIL interferences on three hemostasis analyzers and contribution of a preanalytical module for routine coagulation assays. *Clin Lab*. 2016;62:1979-1987.
 68. Dimeski G, Jones BW. Lipemic samples: effective process for lipid reduction using high speed centrifugation compared with ultracentrifugation. *Biochem Med*. 2011;21:86-94.
 69. Negrini D, Bernadi D, Antonelli G, Plebani M. Interference of lipemia in samples for routine coagulation testing using a Sysmex CS-5100 coagulometer. *Int J Lab Hematol*. 2019;41:772-777.
 70. Gardiner C, Lane P, Taylor H, Mackie IJ. A practical method for reducing the interference due to lipemia in coagulation tests. *Int J Lab Hematol*. 2020;42:140-144.
 71. van Rossum AP, Vlasveld LT, van den Hoven LJ, de Wit AWM, Castel A. False prolongation of the APTT in inflammatory patients: interference of C-reactive protein. *Br J Haematol*. 2012;57:381-440.
 72. Schouwers SM, Delanghe JR, Devreese KM. Lupus Anticoagulant (LAC) testing in patients with inflammatory status: does C-reactive protein interfere with LAC test results? *Thromb Res*. 2010;125:102-104.
 73. Lippi G, Salvagno GL, Adcock DM, Gelati M, Guidi GC, Favaloro EJ. Right or wrong sample received for coagulation testing? Tentative algorithms for detection of an incorrect type of sample. *Int J Lab Hematol*. 2010;32:132-137.
 74. Gosselin RC, Honeychurch K, Kang HJ, Dwyre DM. Effect of multiple freeze-thaw cycles on coagulation testing. *Semin Thromb Hemost*. 2020;46:515-520.
 75. Woodhams B, Giradot O, Blanco M, Collesse G, Gourmelin Y. Stability of coagulation proteins in frozen plasma. *Blood Coagul Fibrinolysis*. 2001;12:229-236.
 76. Gosselin RC, Dwyre DW. Determining the effect of freezing on coagulation testing: comparison of results between fresh and once frozen-thawed plasma. *Blood Coagul Fibrinolysis*. 2015;26:69-74.
 77. Vis R, van den Besselaar AMHP, Witteveen E, Meeuwse-Braun J, Schafer-van MH. Artificial prolongation of the prothrombin time of lyophilised plasma induced by transportation with solid carbon dioxide. *Thromb Haemost*. 1992;1:725.
 78. Trondsetås L, Mikkelsen G, Lian IA. The effects of dry ice exposure on plasma pH and coagulation analyses. *Clin Chem Lab Med*. 2017;56:59-64.
 79. Odsæter IH, Lian IA, Bratberg K, Mikkelsen G. Dry ice exposure of plasma samples influences pH and lupus anticoagulant analysis. *Chem Lab Med*. 2015;53:809-813.

80. Jo S-Y, Kim I-S, Kim H. Potential misdiagnosis of von Willebrand disease caused by insufficient mixing after thawing of frozen plasma. *Int J Lab Hematol*. 2020;43:831-836.
81. Zhao Y, Feng G, Feng L. Effects of pre analytical storage time, temperature and freeze-thaw time on coagulation factor activities in citrate-anticoagulated plasma. *Ann Transl Med*. 2018;6:456-465.
82. Maelegheer K, Devreese KMJ. The impact of repeated freeze-thaw cycles on antiphospholipid antibody titer. *Res Pract Thromb Haemost*. 2018;6:366-369.
83. Clinical and Laboratory Standards Institute. *Platelet Function Testing by Aggregometry: Approved Guideline*. Clinical and Laboratory Standards Institute; 2008;28:1-45.
84. Harrison P, Mackie I, Mumford A, et al. Guidelines for the laboratory investigation of heritable disorders of platelet function. *Br J Haematol*. 2011;155:30-44.
85. Cattaneo M, Cerletti C, Harrison P, et al. Recommendations for the Standardization of Light Transmission Aggregometry: A Consensus of the Working Party from the Platelet Physiology Subcommittee of SSC/ISTH. *J Thromb Haemost*. 2013;11:1183-1189.
86. Le Quellec S, Bordet J-C, Negrier C, Dargaud Y. Comparison of current platelet functional tests for the assessment of aspirin and clopidogrel response. A review of the literature. *Thromb Haemost*. 2016;116:638-650.
87. Kottke-Marchant K, Corcoran G. The laboratory diagnosis of platelet disorders. *Arch Pathol Lab Med*. 2002;126:133-146.
88. Germanovich K, Femia EA, Cheng CY, Dovlatova N, Marco CM. Effects of pH and concentration of sodium citrate anticoagulant on platelet aggregation measured by light transmission aggregometry induced by adenosine diphosphate. *Platelets*. 2018;29:21-26.
89. Stegnar M, Knezevic A, Bozic-Mijovski M. The effect of pre-analytical variables on light transmittance aggregometry in citrated platelet-rich plasma from healthy subjects. *Clin Chem Lab Med*. 2010;48:1463-1465.
90. Andreasen JB, Pistor-Riebold TU, Knudsen IH, Ravn HB, Anne-Hvas AM. Evaluation of different sized blood sampling tubes for thromboelastometry, platelet function, and platelet count. *Clin Chem Lab Med*. 2014;52:701-706.
91. Hellstern P, Sturzebecher U, Wuchold B, et al. Preservation of in vitro function of platelets stored in the presence of a synthetic dual inhibitor of factor Xa and thrombin. *J Thromb Haemost*. 2007;5:2119-2126.
92. Haubelt H, Vogt A, Hellstern P. Preservation of platelet aggregation and dense granule secretion during extended storage of blood samples in the presence of a synthetic dual inhibitor of factor Xa and thrombin. *Platelets*. 2008;19:496-501.
93. Tauer JT, Lohse J, König S, Knöfler R. Investigations of platelet function in whole blood with BAPA as anticoagulant. *Transfus Med Hemother*. 2010;37:284-288.
94. Mani H, Hellis M, Lindhoff-Last E. Platelet function testing in hirudin and BAPA anticoagulated blood. *Clin Chem Lab Med*. 2011;49:501-507.
95. Chapman K, Favaloro EJ. Time dependent reduction in platelet aggregation using the multiplate analyser and hirudin blood due to platelet clumping. *Platelets*. 2018;29:305-308.
96. Favaloro EJ. Utility of the PFA-100 for assessing bleeding disorders and monitoring therapy: a review of analytical variables, benefits and limitations. *Haemophilia*. 2001;7:170-179.
97. Jilma B. Platelet function analyser (PFA-100): a tool to quantify congenital or acquired platelet dysfunction. *J Lab Clin Med*. 2001;138:152-163.
98. Heilmann EJ, Kundu SK, Sio R, Garcia C, Gomez R, Christie DJ. Comparison of four commercial citrate blood collection systems for platelet function analysis by the PFA-100 system. *Thromb Res*. 1997;87:159-164.
99. Harrison P, Robinson MSC, Mackie IJ, et al. Performance of the platelet function analyser PFA-100R in testing abnormalities of primary haemostasis. *Blood Coagul Fibrinolysis*. 1999;10:25-31.
100. Welch EL, Crooks MG, Hart SP. Agreement between blood draw techniques for assessing platelet activation by flow cytometry. *Platelets*. 2019;30:530-534.
101. Schmitz G, Rothe G, Ruf A, et al. European working group on clinical cell analysis: consensus protocol for the flow cytometric characterisation of platelet function. *Thromb Haemost*. 1998;79:885-896.
102. Goodall AH, Appleby J. Flow-cytometric analysis of platelet membrane glycoprotein expression and platelet activation. *Methods Mol Biol*. 2004;272:225-253.
103. Michelson AD, Linden M, Barnard MR, Furman MI, Frelinger AL. Flow cytometry. In: Michelson AD, ed. *Platelets*. Academic Press; 2007: 545-563.
104. Lance MD, Henskens YMC, Nelemans P, et al. Do blood collection methods influence whole-blood platelet function analysis? *Platelets*. 2013;24:275-281.
105. Le Quellec S, Paris M, Nougier C, et al. Pre-analytical effects of pneumatic tube system transport on routine haematology and coagulation tests, global coagulation assays and platelet function assays. *Thromb Res*. 2017;153:7-13.
106. Poletaev AV, Koltsova EM, Ignatova AA, et al. Alterations in the parameters of classic, global, and innovative assays of hemostasis caused by sample transportation via pneumatic tube system. *Thromb Res*. 2018;170:156-164.
107. Wallin O, Soderberg J, Grankvist K, Andreas Jonsson P, Hultdin J. Pre-analytical effects of pneumatic tube transport on routine haematology, coagulation parameters, platelet function and global coagulation. *Clin Chem Lab Med*. 2008;46:1443-1449.
108. Thalen S, Forsling I, Eintrei J, Soderblom L, Antovic JP. Pneumatic tube transport affects platelet function measured by multiplate electrode aggregometry. *Thromb Res*. 2013;132:77-80.
109. Femia EA, Pugliano M, Podda G, Cattaneo M. Comparison of different procedures to prepare platelet-rich plasma for studies of platelet aggregation by light transmission aggregometry. *Platelets*. 2012;23:7-10.
110. Merolla M, Nardi MA, Berger JS. Centrifugation speed affects light transmission aggregometry. *Int J Lab Hematol*. 2012;34:81-85.
111. Cattaneo M, Lecchi A, Zighetti ML, Lussana F. Platelet aggregation studies: autologous platelet-poor plasma inhibits platelet aggregation when added to platelet-rich plasma to normalize platelet count. *Haematologica*. 2007;92:694-697.
112. Cadamuro J, Simundic AM, Ajner E, Sandberg S. A pragmatic approach to sample acceptance and rejection. *Clin Biochem*. 2017;50:579-581.

How to cite this article: Kitchen S, Adcock DM, Dauer R, et al. International Council for Standardization in Haematology (ICSH) recommendations for processing of blood samples for coagulation testing. *Int J Lab Hematol*. 2021;43:1272-1283. <https://doi.org/10.1111/ijlh.13702>

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